Instructions for Use

Caliper Rings

Morcher GmbH

Kapuzinerweg 12 | 70374 Stuttgart | Germany

www.morcher.com

CE 0483

	[en] ENGLISH INSTRUCTIONS FOR USE CALIPER RINGS.....	3
	[de] DEUTSCH GEBRAUCHSANWEISUNG CALIPER RINGE	11
	[es] ESPAÑOL INSTRUCCIONES DE USO ANILLOS DEMARCADORES	20
	[et] EESTI KEEL KASUTUSJUHEND NIHIKURÕNGAD	29
	[fr] FRANÇAIS NOTICE ANNEAUX DE CALIBRAGE	37
	[hr] HRVATSKI UPUTE ZA UPORABU PRSTENASTE MJERKE.....	46
	[it] ITALIANO ISTRUZIONI PER L'USO ANELLI DI CALIBRAZIONE	54
	[lt] LIETUVIŲ NAUDOJIMO INSTRUKCIJOS PALAIKIMO ŽIEDAI	62
	[nl] NEDERLANDS GEBRUIKSAANWIJZINGEN MEETRINGEN	70
	[pl] POLSKI INSTRUKCJA UŻYTKOWANIA PIERŚCIEŃ KALIBRACYJNE CALIPER RING	79
	[ro] ROMÂNĂ INSTRUȚIUNI DE UTILIZARE INELE MARCATOARE	88
	[sv] SVENSKA BRUKSANVISNING MÄTRINGAR	97

MODELS

REF

4, 4L, 5, 6, 76E

DISCLAIMER

Please note: Not all models are available in all countries and from all distributors.

PRODUCT DESCRIPTION

The Caliper Ring is a sterile medical device to enable the ophthalmologists/ophthalmic surgeons to create a more precise anterior continuous curvilinear capsulorhexis (ACCC) during intraocular lens (IOL) implantation. It is a closed ring, made from one piece of black polyethylmethacrylate (PEMA Black) and available in different sizes.

The Caliper Ring is an accessory designed to size and centre the position of the anterior capsulorhexis. The rings are available in different sizes. The Caliper Ring should be removed after the anterior continuous curvilinear capsulorhexis (ACCC).

SPECIFICATIONS

Model (REF)	Total Ø [mm]	ACCC [mm]
4	4.3	4.0
4L	4.8	4.5
5	5.3	5.0
6	6.3	6.0
76E	6.8	6.5

ACCC = anterior continuous curvilinear capsulorhexis

INTENDED PURPOSE

The Caliper Ring is intended as accessory for the implantation of a Bag-In-the-Lens (BIL) and BioComFold IOL (BCF). The Caliper Ring is designed as aid for the surgical step of performing an ACCC.

The Caliper Ring is available in different sizes corresponding to the required capsulorhexis size.

INTENDED USER

The Caliper Ring is used by either an ophthalmic surgeon or an ophthalmologist during a state-of-the-art operation.

USE ENVIRONMENT

The Caliper Ring is used under aseptic conditions in an operating room or a doctor's practice using a surgical microscope.

TARGET GROUP

Children and adults of any age (possibly limited by the target group of the IOL).

The Caliper Ring has not been clinically investigated in pregnant or nursing women. Ophthalmologists/ophthalmic surgeons, considering usage in such patients, should carefully assess the benefit-risk ratio before usage.

INDICATIONS

The use of the Caliper Ring is advisable in any surgical procedure where a specified ACCC is desired.

CONTRAINDICATIONS

When used as intended, there are no contraindications for the use of the Caliper Ring.

INTENDED CLINICAL BENEFIT

The Caliper Ring can be used as an accessory for the implantation of BCF IOL, including BIL IOL. It provides the benefit of a precisely sized and positioned anterior continuous curvilinear capsulorhexis, facilitating safe and effective implantation of the IOL.

PERFORMANCE CHARACTERISTICS

The primary performance parameter of the Caliper Rings is technical success, which is implicitly shown in the scientific literature by success of the IOL implantation procedure. Due to the lack of reported Caliper Ring-related problems and complications, a technical success rate of 100 % can be assumed.

LIFETIME

The Caliper Ring is removed after the rhexis.

SIDE EFFECTS

At the present time, there are no known side effects.

WARNINGS

- Do not use if the sterile packaging is already open or damaged.
- Do not re-sterilise the Caliper Ring. Re-sterilization can lead to damage to the material and loss of functionality of the device.
- Do not re-use the Caliper Ring. The Caliper Ring is intended for single use. Re-use can lead to reduced effectiveness or functionality of the Caliper Ring and/or loss of sight.
- Do not use the Caliper Ring if the expiration date has expired.
- Rinse the Caliper Ring only with sterile rinsing solutions such as sterile Ringer's solution or sterile BSS® (balanced salt solution).
- Do not use a damaged device. A Caliper Ring that is damaged during usage must be disposed of immediately.
- Improper handling and/or folding may lead to damage of the Caliper Ring and may cause injury.
- Observe accompanying documents and warning instructions before using the device.
- Consult Instructions for Use before using the device.

PRECAUTIONARY MEASURES

QUALIFICATION AND TRAINING FOR USERS

- The Caliper Ring may only be used by an ophthalmologist/ophthalmic surgeon.
- The device may only be used by a doctor experienced in the implantation of IOL, or under his/her supervision. The use of the device requires a detailed understanding of the technical fundamental principles, clinical applications and the risks of the device.
- Before using the Caliper Ring, the surgeon should carefully read the Instructions for Use and all the information material provided by the manufacturer for the proper handling and use of the device.

STERILISATION AND PACKAGING

- The Caliper Ring is sterilised by gamma irradiation.
- The Caliper Ring is supplied dry and sterile.
- The primary container is sealed, germproof, and packed into a sealed and germproof transparent pouch (double sterile barrier system).
- Sterility is only guaranteed as long as the transparent pouch and the primary container are neither opened nor damaged.
- The expiration date is printed on the folded box and the primary container.

STORAGE CONDITIONS / TRANSPORT CONDITIONS

- Store flat at room temperature + 15 °C up to + 25 °C.
- Transport at - 18 °C up to + 60 °C.
- Protect from direct sunlight and extreme temperatures.
- Protect from moisture/water and store in a dry environment.

PATIENT INFORMATION

- Each patient shall be fully informed about the device, especially about potential risks (including any undesirable side effects, warnings, precautions, and contraindications) and benefits, before implantation.
- The surgeon performing the surgery shall inform the patient of all side effects and risks, particularly those that have not already been pointed out in the initial consultation.
- The patient shall be instructed to inform the doctor about all side effects, which he/she may observe after the surgery.
- The Caliper Ring is removed after the rhexis and is therefore neither exposed to specific environmental conditions nor medical procedures or examinations like laser treatment, magnetic resonance imaging (MRI) or X ray scan. Regardless, such conditions or procedures and examinations would not negatively impact the performance of the device and would not pose any risk to the device.
- The patient shall be informed that there are no known circumstances, in which the device interfere with certain diagnostic tests, evaluations, or therapeutic treatments or other procedures, or interfere with the equipment used in these procedures.

NOTES ON HANDLING

1. Check all labels with regard to type, specific details and expiration date before opening the transparent pouch.
2. Only pick up the Caliper Ring with smooth-tipped forceps and handle it carefully to avoid damage.
3. Under aseptic conditions, remove the device as follows:
 - Open the transparent pouch at the points indicated and remove the primary container.
 - Carefully pull off the seal from the primary container.
 - Carefully remove the Caliper Ring; do not pull on the Caliper Ring.
4. Check Caliper Ring before use and only use an undamaged Caliper Ring.
5. It is recommended to rinse the Caliper Ring with sterile Ringer's solution or sterile BSS® prior to usage.

6. PERFORMING THE ACCC

The anterior chamber must be properly filled with viscoelastic before introducing the Caliper Ring. The Caliper Ring can be grasped with two forceps at the 3 and 9 o'clock position and formed into an oval shape. This device can then easily be inserted through a 3.0 mm incision into the anterior chamber, which has been filled with viscoelastic. The Caliper Ring is positioned on the anterior capsule by using a hook. In order to stabilise the Caliper Ring, once positioned on top of the anterior capsule, more viscoelastic can be injected until the anterior surface of the lens capsule expresses a meniscus, which can have the same diameter as the Caliper Ring. The ACCC is initiated in the centre of the Caliper Ring and performed by tearing the anterior capsule along the inner border of the Caliper Ring. Once the ACCC is completed, the Caliper Ring is removed, e.g. using capsulorhexis forceps.

USE WITH OTHER DEVICES AND / OR GENERAL PURPOSE EQUIPMENT

The Caliper Ring can be used as an accessory for MORCHER® BioComFold IOL and BIL IOL.

There are no restrictions known regarding the combination of the device with general surgical equipment, ophthalmic viscosurgical devices (OVD, viscoelastics), or intraocular physiological rinsing solutions.

DISPOSAL

- The packaging and its contents shall be disposed of in accordance with the current local regulations.
- There are no special disposal regulations applicable for devices, which are uncontaminated.
- Contaminated devices/packaging shall be handled and disposed of in accordance with the provisions of the respective clinic.

EXCLUSION OF LIABILITY

Morcher GmbH and its distributors cannot be held liable for any procedures, methods, or surgical techniques used by the surgeon, nor for the selection of the device for any eye or surgical procedure. This is solely the responsibility of the surgeon.

COMPLAINTS

Complaints should be addressed to

Morcher GmbH

Quality Assurance Department
Kapuzinerweg 12
70374 Stuttgart
Germany

Phone + 49 (0) 711/95 320 - 0
Fax + 49 (0) 711/95 320 - 80
E-Mail complaints@morcher.com

or to your distributor.

REPORTING ON ADVERSE REACTIONS

Any serious incident that occurred in relation to the device should be reported to

Morcher GmbH

Quality Assurance Department
Kapuzinerweg 12
70374 Stuttgart
Germany

Phone + 49 (0) 711/95 320 - 0
Fax + 49 (0) 711/95 320 - 80
E-Mail complaints@morcher.com

and to the competent authority of the member state in which the user and/or patient is established.
Please also inform your distributor if necessary.

RETURN OF A DAMAGED / NON-STERILE DEVICE

If a Caliper Ring is returned, the serial number, batch code and reason for its return must be provided. The packaging, containing the returned item, must be in accordance to the condition of the device and, where necessary, must warn of any contamination.

ABBREVIATIONS AND PACKAGING SYMBOLS

ABBREVIATION / SYMBOL	DESCRIPTION
IOL	Intraocular Lens
CTR	Capsular Tension Ring
PC	Posterior Chamber
AC	Anterior Chamber
ACCC	Anterior Continuous Curvilinear Capsulorhexis
PCCC	Posterior Continuous Curvilinear Capsulorhexis
ICCE	Intracapsular Cataract Extraction
ECCE	Extracapsular Cataract Extraction
PCO	Posterior Chamber Opacification
D	Diopter in aqueous humour
ØB	Body diameter
ØT	Total diameter
ØO	Optic diameter
ØI	Inner diameter
R	Radius
Th	Haptic side view
A-CON	A-constant
AC DEPTH / ACD	Anterior Chamber Depth
PMMA	Polymethylmethacrylate
PEMA	Polyethylmethacrylate
UV	Ultraviolet
HEV	High Energy Visible light
C-F-M	Compression forged method
CQ	Clinical Quality
BSS®	Balanced Salt Solution
SN	Serial Number

ABBREVIATION / SYMBOL	DESCRIPTION
	Unique Device Identification
GTIN	Global Trade Item Number
	Medical Device
	Catalogue number
	Batch code
	Use by: Year - Month
EXP	Use by: Year - Month
	Manufacturer
DE	Country and date of manufacture
	Importer
	Distributor
XXXX	CE Mark with 4-digit number of the notified body
	Sterilised using irradiations
	Sterilised using steam
	Double sterile barrier system
	Single sterile barrier system
kGy	Kilogray
	Do not reuse
	Do not resterilize
	Consult instructions for use
 ifu.morcher.com	Consult electronic instructions for use
	Caution, consult accompanying documents
	Keep dry
	Keep away from sunlight
	Temperature limit (The applicable lower and upper temperature limits for storage are indicated on the device packaging. The applicable lower and upper temperature limits for transport are indicated on the outer packaging used for transportation.)

ABBREVIATION / SYMBOL	DESCRIPTION
	Do not use if packaging is damaged
	This side up
	Recyclable packaging
	Date of implantation
	Healthcare institution
	Patient information website
	Patient identification

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