Instructions for Use

Bag-In-the-Lens

Foldable Intraocular Lenses

Morcher GmbH

Kapuzinerweg 12 | 70374 Stuttgart | Germany

www.morcher.com

CE 0483

	[en] ENGLISH INSTRUCTIONS FOR USE	
	Bag-In-the-Lens Foldable Intraocular Lenses	3
	[de] DEUTSCH GEBRAUCHSANWEISUNG	
	Bag-In-the-Lens Faltbare Intraokularlinsen	15
	[es] ESPAÑOL INSTRUCCIONES DE USO	
	Bag-In-the-Lens Lentes intraoculares plegables	27
	[et] EESTI KEEL KASUTUSJUHEND	
	Bag-In-the-Lens volditavad intraokulaarsed läätsed.....	39
	[fr] FRANÇAIS NOTICE	
	Bag-In-the-Lens Lentille intraoculaire souple.....	51
	[hr] HRVATSKI UPUTE ZA UPORABU	
	Bag-In-the-Lens Savitljive intraokularne leće	63
	[it] ITALIANO ISTRUZIONI PER L'USO	
	Bag-In-the-Lens Lenti intraoculari pieghevoli	74
	[lt] LIETUVIŲ NAUDOJIMO INSTRUKCIJOS	
	Bag-In-the-Lens Sulankstomi intraokuliniai lęšiai	86
	[nl] NEDERLANDS GEBRUIKSAANWIJZING	
	Bag-In-the-Lens vouwbare intraoculaire lenzen.....	97
	[pl] POLSKI INSTRUKCJA UŻYTKOWANIA	
	Bag-In-the-Lens Zwijalne soczewki wewnętrzgałkowe.....	109
	[ro] ROMÂNĂ INSTRUȚIUNI DE UTILIZARE	
	Bag-In-the-Lens Cristaline intraoculare pliabile	121
	[sv] SVENSKA BRUKSANVISNING	
	Bag-In-the-Lens vikbara intraokulära linser	133

MODELS

REF

89A, 89D, 89F

DISCLAIMER

Please note: Not all models are available in all countries and from all distributors.

PRODUCT DESCRIPTION

The foldable Bag-In-the-Lens (BIL) intraocular lens (IOL), hereafter referred to as lens or implant, is a sterile optical medical implant designed to replace the natural lens. The quality of the material used permits folding the optic for small incision insertion. The lens is made from one piece of a hydrophilic acrylic copolymer including a copolymerized UV absorber. The lens consists of a central round optic surrounded by a peripheral groove defined by the elliptical haptics, perpendicularly oriented to each other in order to accommodate the anterior and posterior capsulorhexis.

The lens is available with different sizes of optics in order to ensure optimal fitting of the implant in the eye. Each lens is tested separately for its optical quality and refractive power.

DEVICE MATERIALS

The BIL IOL are made from a hydrophilic acrylic copolymer. This polymer contains < 1 % of a copolymerizable UV absorber.

SPECIFICATIONS

Model (REF)	Total Ø [mm]	Optic Ø [mm]	Recommended Caliper Ring for ACCC
89A	7.5	5.0	model 5, model 4L (children)
89D	6.5	4.5	model 4L
89F	8.5	5.0	model 5

ACCC = anterior continuous curvilinear capsulorhexis

Refractive index = 1.46 for all lenses.

INTENDED PURPOSE

MORCHER® foldable Bag-In-the-Lens (BIL) is intended to serve as replacement of the natural lens and is designed for fixation in the capsular bag. It is used after extra-capsular cataract extraction.

INTENDED USER

The BIL is used by either an ophthalmic surgeon or an ophthalmologist trained on the BIL.

USE ENVIRONMENT

The BIL is used under aseptic conditions in an operating room or a doctor's practice using a surgical microscope.

TARGET GROUP

Target group are children and adults of any age.

The implants have not been clinically investigated in pregnant or nursing women. Ophthalmologists/ ophthalmic surgeons, considering implantation in such patients, should carefully assess the benefit-risk ratio before implantation.

INDICATIONS

- For senile cataract
- As correction of an aphakic eye where other optical aids are not available, or are not implementable
- For traumatic cataract
- In congenital or childhood-acquired cataract

CONTRAINDICATIONS

The usual contraindications for a cataract operation apply here, as well as the following ones:

- Pre-operative complications of cataract surgery (vitreous prolapse, hemorrhage)
- Anticipated post-operative complications of cataract surgery
- Previous intraocular operation on the eye scheduled for surgery
- Combined surgical intervention at the time of implantation, e.g. glaucoma surgery or keratoplasty (iridectomy is acceptable)
- Acute infection or inflammation of the eye
- Chronic use of steroids, immunosuppressant and/or antineoplastic agents
- Significant corneal disease or surface pathology of the eye including sicca syndrome, keratoconus or blepharitis
- Chronic uveitis, iritis, iridocyclitis or rubeosis iridis
- Iris atrophy
- Aniridia
- Proliferative diabetic retinopathy
- Rubeola cataracts
- Corneal dystrophy, particularly endothelial dystrophy
- Severe microphthalmia or macrophthalmia
- Decompensated chronic glaucoma
- Systemic diseases with ophthalmological involvement (e.g. diabetes, complications of AIDS) that might cause misinterpretation of postoperative findings
- Patients with a history of recurrent anterior or posterior segment inflammations of unknown etiology
- Patients with a shallow anterior segment, e.g. microphthalmia or certain types of chronic angle closure glaucoma
- Patients in whom an intraocular lens might interfere with the ability to evaluate, diagnose or treat posterior segment diseases
- Implantation in the functionally better eye, or in patients having only one eye

Certain operative complications may lead to a contraindication against the use of the device, such as:

- Rupture of posterior capsular bag
- Detachment of Descemet's membrane
- Significant anterior chamber bleeding
- Persistent bleeding
- Iris damage
- Vitreous loss
- Irremovable vitreous bulge in the anterior chamber
- Uncontrollable intraocular pressure

INTENDED CLINICAL BENEFIT

The BIL is used to replace the natural lens in case of cataract or aphakia. It leads to the benefit of a restored or improved visual acuity. Due to the material properties, the lens can be folded and implanted through a small incision.

PERFORMANCE CHARACTERISTICS

The BIL IOL significantly improves visual acuity. It fulfils the requirements of the DIN EN ISO 11979-7 on clinical investigations on IOLs: In a clinical study by Tassignon et al. a postoperative decimal CDVA (corrected distance visual acuity) of 0.5 or better was achieved in 99.36 % of eyes without ocular comorbidity and in 89.08 % of all eyes [Tassignon MJ et al., JCRS 2011, 37(12):2120]. For details see section "Summary of safety and clinical performance".

LIFETIME

The lifetime of the BIL is 10 years. Without any complications or adverse events, the device may stay in the patient for his lifetime.

COMPLICATIONS OF THE IMPLANTATION (SPECIFIC AND NON-SPECIFIC)

- Complications associated with the general surgical procedure
- Shallow anterior chamber
- Corneal oedema
- Corneal dystrophy
- Pupillary block
- Intrapupillary membrane
- Iridocyclitis, vitritis, uveitis
- Precipitates on the intraocular lens (calcification/opacification of the intraocular lens)
- Decentration of the intraocular lens
- (Sub-) luxation of the intraocular lens
- Dislocalization of the intraocular lens
- Iris atrophy
- Infection, endophthalmitis and panophthalmitis
- Suprachoroidal hemorrhage
- Hyphema
- Anterior synechia
- Uncontrolled elevated intraocular pressure

- Glaucoma
- Vitreous bulge
- Vitreous detachment
- Retinal detachment
- Visual axis opacification (VAO)/post-cataract membrane/secondary membrane
- Inflammatory pupillary membrane
- Cystoid macular oedema
- Anterior capsular contraction syndrome (ACCS)
- Capsular block syndrome
- Toxic anterior segment syndrome (TASS)
- Hypopyon
- Undesired ametropia
- Aniseikonia
- Surgically induced astigmatism
- Iris capture
- Dysphotopsia; glare

SIDE EFFECTS

There have been occasional reports of secondary glaucoma after the implantation of posterior chambers lenses in patients with controlled glaucoma. (In pediatric cataract surgery, especially in patients aged ≤ 12 weeks at surgery, higher incidences of secondary glaucoma occur, which are related to the procedure of lensectomy in general.)

Occasionally, as a result of intra-operative complications of the cataract extraction, post-operative problems such as macular edema, hyphema or vitritis occur, so such patients should be carefully monitored post-operatively.

WARNINGS

- Do not use if the sterile packaging is already open or damaged.
- Do not re-sterilise the implant. Re-sterilization can lead to damage to the material and loss of functionality of the implant.
- Do not re-use the implant. The implant is intended for single use. Re-use can lead to reduced effectiveness or functionality of the implant and/or loss of sight.
- Do not use the implant if the expiration date has expired.
- Rinse the implant only with sterile rinsing solutions such as sterile Ringer's solution or sterile BSS® (balanced salt solution).
- Do not use a damaged implant. An implant that is damaged during the implantation must be disposed of immediately.
- Improper handling and/or folding of the implant may lead to damage of the implant and may cause injury.
- Observe accompanying documents and warning instructions before using the device.
- Consult Instructions for Use before using the device.
- Do not use implants from dried-out packages.
- The implant must remain moistened while exposed to the air.
- Do not use the implant if it may have been contaminated with unsterile material (e.g. if it falls on the floor). Dispose of the potentially contaminated implant and use a new product.

PRECAUTIONARY MEASURES

QUALIFICATION AND TRAINING FOR USERS

- The lens may only be used by an ophthalmologist/ophthalmic surgeon trained on the BIL as specified by Prof. Dr. em. Tassignon. Please contact Morcher GmbH for more information.
- The use of this product requires a thorough understanding of the technical fundamental principles, clinical application and the risks of using this product.
- The surgeon performing the implantation should carefully read all the information material provided by the manufacturer for the proper use and insertion of the implant before carrying out the implantation.

CALCIFICATION OF THE IMPLANT

In rare cases opacification/calcification of the IOL may occur, which usually results in loss of vision. The exact pathomechanism has not yet been fully understood. To minimize the risk of IOL calcification, the IOL should not come into contact with phosphate-containing substances. We thus recommend avoiding phosphate-containing substances during surgery.

NOTE WHEN USING LASERS

The laser beam must be focused directly on the target point behind the implant. If a laser beam is focused on the implant, it can result in damage to the implant.

LIMITATIONS

A positive outcome of IOL implantation cannot be assured in patients with pre-existing ocular disorders (chronic medicinal miosis, diabetic retinopathy, previous cornea transplant, history of retinal detachment or iritis, glaucoma, corneal disease etc.). Ophthalmologists/ophthalmic surgeons, considering lens implants in such patients, should explore the use of alternative methods. Close postoperative monitoring of intraocular pressure is indicated in patients with pre-existing ocular disorder.

It should be kept in mind that central visual acuity may not be improved in patients with age-related macular or retinal degeneration.

An IOL with a small optic may result in a patient experiencing glare under certain light conditions with large pupils. Ophthalmologists/ophthalmic surgeons should consider this before implanting an IOL having a small optic.

In the case of patients with an impaired calcium-phosphate balance (e.g. patients with diabetes), there is an increased risk of the IOL possibly becoming more cloudy at a later point (calcification).

Colouring substances in connection with hydrophilic implants should be avoided whereas this may cause permanent colouring of the implant material.

STERILISATION AND PACKAGING

- The implant is sterilised by pressurized steam (autoclave).
- The implant is supplied wet and sterile.
- The primary container is sealed, germproof, and packed into a sealed and germproof transparent pouch (double sterile barrier system).
- Sterility is only guaranteed as long as the transparent pouch and the primary container are neither opened nor damaged.
- The expiration date is printed on the folded box and the primary container.

STORAGE CONDITIONS / TRANSPORT CONDITIONS

- Store flat at room temperature + 15 °C up to + 25 °C.
- Transport at - 18 °C up to + 60 °C.
- Protect from direct sunlight and extreme temperatures.
- Protect from moisture/water and store in a dry environment.

PATIENT INFORMATION

- Each patient shall be fully informed about the implant, especially about potential risks (including any undesirable side effects, warnings, precautions, contraindications, and limitations) and benefits, before implantation.
- The surgeon performing the implantation shall inform the patient of all side effects and risks, particularly those that have not already been pointed out in the initial consultation.
- The patient shall be instructed to inform the doctor about all side effects, which he/she may observe after the implantation, especially in case of malfunction of the device or changes in its performance.
- The patient shall be informed that there are no specific environmental conditions known that can affect the BIL IOL. Medical procedures or examinations like laser treatment, magnetic resonance imaging (MRI) or X ray scan do not negatively impact the performance of the device and do not pose any risk to the device.
- The patient shall be informed that there are no known circumstances, in which the device interfere with certain diagnostic tests, evaluations, or therapeutic treatments or other procedures, or interfere with the equipment used in these procedures. Please note, that implantation is contraindicated in patients, in whom the IOL might interfere with the ability to evaluate, diagnose or treat posterior segment diseases.
- Written information about the devices intended for the patient is available at <https://patient-information.morcher.com>.

IMPLANT CARD

The relevant details and one label shall be attached onto the enclosed implant card. The date of the implantation shall be written on this card and the eye that was operated on shall be crossed off. Instructions for filling out the implant card are available at <https://ifu.morcher.com>.

The completed card shall be given to the patient for safekeeping to ensure traceability of the implant; it should be presented to every eye specialist whom the patient consults

NOTES ON HANDLING / IMPLANTATION

1. Check all labels with regard to type, specific details and expiration date before opening the transparent pouch.
2. Only pick up the implant with smooth-tipped forceps and handle it carefully to avoid damage.
3. Under aseptic conditions, remove the implant as follows:
 - Open the transparent pouch at the points indicated and remove the primary container
 - Carefully pull off the seal from the primary container
 - Carefully remove the implant. Do not pull on the implant.
4. Check implant before use and only use an undamaged implant.
5. It is recommended to rinse the implant with sterile Ringer's solution or sterile BSS® prior to the implantation.

6. Folding

- It is important to coat the lens with sufficient viscoelastic before folding with forceps.
- Lens design and power must be considered in the selection of forceps, and the incision size must be selected accordingly.

After implantation of the lens, please ensure that the lens is completely unfolded and correctly positioned.

A secondary iridectomy for pupillary block may be avoided if one or more iridectomies are performed at the time of intraocular lens implantation.

USE WITH OTHER DEVICES AND / OR GENERAL PURPOSE EQUIPMENT

The lenses can be used with the accessory MORCHER® Caliper Ring. The Caliper Ring is an accessory designed to size and centre the position of the anterior capsulorhexis. Caliper Rings are available in different sizes. The size of the Caliper Ring has to match with the groove of the lens (details can be found in the Table above).

There are no restrictions known regarding the combination of the lenses with general surgical equipment, ophthalmic viscosurgical devices (OVD, viscoelastics), or intraocular physiological rinsing solutions.

To minimize the risk of IOL calcification, phosphate-containing substances should be avoided during surgery.

DISPOSAL

- The packaging and its contents shall be disposed of in accordance with the current local regulations.
- There are no special disposal regulations applicable for devices, which are uncontaminated.
- Contaminated devices/packaging shall be handled and disposed of in accordance with the provisions of the respective clinic.

EXCLUSION OF LIABILITY

Morcher GmbH and its distributors cannot be held liable for any procedures, methods, or surgical techniques used by the surgeon, nor for the selection of the implant for any eye. This is solely the responsibility of the surgeon.

COMPLAINTS

Complaints should be addressed to

Morcher GmbH

Quality Assurance Department
Kapuzinerweg 12
70374 Stuttgart
Germany

Phone + 49 (0) 711/95 320 - 0
Fax + 49 (0) 711/95 320 - 80
E-Mail complaints@morcher.com

or to your distributor.

REPORTING ON ADVERSE REACTIONS

Any serious incident that occurred in relation to the device should be reported to

Morcher GmbH

Quality Assurance Department
Kapuzinerweg 12
70374 Stuttgart
Germany

Phone + 49 (0) 711/95 320 - 0
Fax + 49 (0) 711/95 320 - 80
E-Mail complaints@morcher.com

and to the competent authority of the member state in which the user and/or patient is established. Please also inform your distributor if necessary.

RETURN OF A DAMAGED / NON-STERILE DEVICE

If an implant is returned, the serial number, batch code and reason for its return must be provided. The packaging, containing the returned item, must be in accordance to the condition of the implant and, where necessary, must warn of any contamination.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The summary of safety and clinical performance (SSCP) referred to in Article 32 of the Medical Device Regulation (EU) 2017/745 is linked to the Basic UDI-DI number 4048509IOL01S3 and is available at <https://ifu.morcher.com> or in the European database on medical devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed> (as soon as this EUDAMED functionality is available).

ABBREVIATIONS AND PACKAGING SYMBOLS

ABBREVIATION / SYMBOL	DESCRIPTION
IOL	Intraocular lens
CTR	Capsular Tension Ring
PC	Posterior Chamber
AC	Anterior Chamber
ACCC	Anterior Continuous Curvilinear Capsulorhexis
PCCC	Posterior Continuous Curvilinear Capsulorhexis
ICCE	Intracapsular Cataract Extraction
ECCE	Extracapsular Cataract Extraction
PCO	Posterior Chamber Opacification
D	Diopter in aqueous humour
ØB	Body diameter
ØT	Total diameter
ØO	Optic diameter
ØI	Inner diameter
R	Radius
Th	Haptic side view
A-CON	A-constant
AC DEPTH/ACD	Anterior Chamber Depth
PMMA	Polymethylmethacrylate
PEMA	Polyethylmethacrylate
UV	Ultraviolet
HEV	High Energy Visible light
C-F-M	Compression forged method
CQ	Clinical Quality
BSS®	Balanced Salt Solution
SN	Serial Number

ABBREVIATION / SYMBOL	DESCRIPTION
UDI	Unique Device Identification
GTIN	Global Trade Item Number
MD	Medical Device
REF	Catalogue number
LOT	Batch code
	Use by: Year - Month
EXP	Use by: Year - Month
	Manufacturer
DE	Country and date of manufacture
	Importer
	Distributor
CE XXXX	CE Mark with 4-digit number of the notified body
STERILE R	Sterilised using irradiations
STERILE	Sterilised using steam
	Double sterile barrier system
	Single sterile barrier system
kGy	Kilogray
	Do not reuse
	Do not resterilize
	Consult instructions for use
ifu.morcher.com	Consult electronic instructions for use
	Caution, consult accompanying documents
	Keep dry
	Keep away from sunlight
	Temperature limit (The applicable lower and upper temperature limits for storage are indicated on the device packaging. The applicable lower and upper temperature limits for transport are indicated on the outer packaging used for transportation.)

ABBREVIATION / SYMBOL	DESCRIPTION
	Do not use if packaging is damaged
	This side up
	Recyclable packaging
	Date of implantation
	Health care centre or doctor
	Patient information website
	Patient identification

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