

PATIENT CHAIR "MALIN"

EXAMINATION



MADE IN SWEDEN

PATIENT CHAIR "MALIN"

- **Easy rolling chairs that can handle heavy loads and have many flexible accessories**
- **Central brake with extra large rubber wheels that makes the chair very stable in brake mode**
- **Electrical control with a wide height adjustment to fit various types of examinations and treatments**
- **Stable design that provides safety for patient as well as personnel**

Rinis battery powered patient chairs in the series "Malin" are ergonomically designed and can handle heavy patients. The height adjustment between 500-800 mm is controlled by a fixed or wireless hand control.

The chair is equipped with central brake, shaped as a hoop between the rear wheels. The hoop breaks both rear wheels simultaneously and makes the chair almost impossible to move when breaks are locked. While as unbraked the chair is easily rolled and comfortable with the four large wheels.

The chairs have accessories such as foot support, driving handle, wireless hand control, phlebotomy armrest.

There is an extra high version of "Malin" with seat height from 540 mm up to 940 mm for ophthalmologists, to optimize the ergonomic work environment during examination and treatment of children.



TECHNICAL INFORMATION

GENERAL

Length (transport mode)	580 mm
Width (transport mode)	530 mm
Weight unit	41 kg/91 lbs
Height adjustment from floor	500-800 mm
Lifting capacity	120 kg/265 lbs

BACK SUPPORT

Length	360 mm
Width	400 mm
Height adjustment from seat	30-100 mm

SEAT

Length	420 mm
Width	460 mm

ELECTRICAL SPECIFICATION

Battery	24 V
Charger	100-240 V



Rini Ergoteknik AB Fågelsångsvägen 3B
SE-186 42 Vallentuna Sweden
Tel +46 8-594 771 70 Fax +46 8-929 279
info@rini.se www.rini.se