

OASIS[®]





EXPAND YOUR VIEW.
SUSTAIN YOUR ACCESS.

OASIS® IRIS EXPANDER

PATENT NO. 8,439,833
PATENT NO. D686,729



A polypropylene ring that gently expands the iris in cases of intraoperative miosis and/or IFIS and sustains visibility throughout the surgical procedure. Each OASIS Iris Expander® is easily inserted and removed. Once removed, the iris returns to its normal shape and function.

- 7.0mm and 6.25mm diameter
- Molded polypropylene ring
- Four pockets where iris is cradled
- Single use, packaged sterile with inserter

Qty	Diameter	Ref #
5 per box	7.0mm	9700
1 per box	7.0mm	9700S
5 per box	6.25mm	9625
1 per box	6.25mm	9625S



SOFT SHIELD® COLLAGEN CORNEAL SHIELD

CE MARK PENDING

COMFORT & PROTECTION

Provides ocular surface protection following surgery, injection and traumatic or non-traumatic corneal conditions. Available in Quick Shield, 12 hour, 24 hour, and 72 hour dissolution times.


Diameter: 14.5mm

Compound base curve: 9mm

- Sterile, 10 per box
- Coats the surface of the eye with collagen film as it gradually dissolves
- Protects the surface of the eye
- Certified suitable for use by EDQM

Description	Ref #
SOFT SHIELD QS Collagen Corneal Shield, Quick Shield	7002
SOFT SHIELD Collagen Corneal Shield, 12 hour	7012
SOFT SHIELD Collagen Corneal Shield, 24 hour	7024
SOFT SHIELD Collagen Corneal Shield, 72 hour	7072

EMPOWER YOUR PATIENTS WITH CHOICE.




AVAILABLE IN THE USA EXCLUSIVELY FROM OASIS® MEDICAL

Intacs corneal implants are indicated for the treatment of keratoconus to eliminate or reduce the myopia and the astigmatism in those patients that are no longer able to achieve satisfactory vision correction with contact lenses or spectacles and want to defer the need for a corneal transplant.

- For the treatment of keratoconus
- FDA approved for the treatment of myopia
- Stabilizes cornea
- Reshapes cornea to more natural state
- Restores contact lens tolerance
- Defers corneal transplants
- Minimally invasive

Verification of IRB Approval required prior to shipment of Intacs Corneal Implants for Keratoconus.

Product Name	Thickness	Symmetry	Ref #
Intacs Corneal Implants for Myopia	0.250mm	Symmetric	11700-250
Intacs Corneal Implants for Myopia	0.300mm	Symmetric	11700-300
Intacs Corneal Implants for Myopia	0.350mm	Symmetric	11700-350
Intacs Corneal Implants for Keratoconus	0.210mm	Symmetric	11600-210
Intacs Corneal Implants for Keratoconus	0.250mm	Symmetric	11600-250
Intacs Corneal Implants for Keratoconus	0.300mm	Symmetric	11600-300
Intacs Corneal Implants for Keratoconus	0.350mm	Symmetric	11600-350
Intacs Corneal Implants for Keratoconus	0.400mm	Symmetric	11600-400
Intacs Corneal Implants for Keratoconus	0.450mm	Symmetric	11600-450
Intacs Corneal Implants for Keratoconus	0.210mm/350mm	Asymmetric	11400-210/350
Intacs Corneal Implants for Keratoconus	0.210mm/400mm	Asymmetric	11400-210/400
Intacs Corneal Implants for Keratoconus	0.210mm/450mm	Asymmetric	11400-210/450





Feather Reusable scalpels are designed for repeated use, can be resterilized at your surgical facility, and are easily disposed of after multiple uses. Feather Disposable and Safeshield scalpels are single-use.

Reusable	■ 5 per box	■ Sure-grip aluminum handle	■ Individual autoclavable cases	
Disposable	■ 5 per box	■ Disposable plastic handle	■ Sterile package	
Safeshield	■ 5 per box	■ Disposable plastic handle	■ Sterile package	■ Retractable safety guard



Controlled Depth Groove

Depth	Reusable
300 microns (0.3mm)	07330G
350 microns (0.35mm)	07335G
400 microns (0.4mm)	07340G
550 microns (0.55mm)	07355G
600 microns (0.6mm)	07360G



Glaucoma Groove

Depth	Thickness	Reusable
400 microns (0.4mm)	0.1mm	06340G



Implant

Width	Bevel	Blade	Reusable
3.5mm	Up	Angled	07635BR
5.1mm	Up	Angled	07651BR



Ultra-Sharp Slit

Width	Bevel	Blade	Reusable	Disposable
1.4mm	Up	Angled	0614BUC	P-0614BUC
1.6mm	Up	Angled	0616BUC	P-0616BUC
1.8mm	Up	Angled	0618BUC	P-0618BUC
2.2mm	Up	Angled	0622BUC	P-0622BUC
2.4mm	Up	Angled	0624BUC	P-0624BUC
2.75mm	Up	Angled	06275BUC	P-0675BUC
2.8mm	Up	Angled	0628BUC	P-0628BUC
3.0mm	Up	Angled	0630BUC	P-0630BUC
3.2mm	Up	Angled	0632BUC	P-0632BUC



Slit

Width	Bevel	Blade	Reusable	Disposable	Safeshield
2.5mm	Double	Angled	07625B		
2.65mm	Double	Angled	07627B		
2.8mm	Double	Angled	07628B	P-7628B	PSS-7628B
3.0mm	Double	Angled	07630B	P-7630B	
3.0mm	Up	Angled	07630BU	P-7630BU	
3.2mm	Double	Angled	07632B	P-7632B	
3.2mm	Up	Angled	07632BU		



Clear Cornea Slit

Width	Bevel	Thickness	Reusable	Disposable
2.8mm	Double	0.15mm	09628B	P-9628B
3.0mm	Double	0.15mm	09630B	P-9630B
3.2mm	Double	0.15mm	09632B	P-9632B



Crescent Tunnel

Width	Bevel	Reusable	Disposable
2.2mm	Double	07922B	
2.6mm	Up	07926BU	P-7926BU
2.6mm	Down	07926BD	
2.6mm	Double		P-7926B



Parallel Slit

Width	Bevel	Thickness	Reusable
2.2mm	Up	0.15mm	07622BUPS



Round Tunnel

Width	Bevel	Blade	Reusable	Disposable	Safeshield
2.2mm	Double	Angled	07722B		
2.6mm	Double	Angled	07726B	P-7726B	
3.0mm	Double	Angled	07730B		PSS-7730B



Stab Incision

Angle	Blade	Reusable	Disposable	Safeshield
15°	Straight	715	P-715	PSS-715
22.5°	Straight	722	P-722	PSS-730
30°	Straight	730	P-730	
45°	Straight	745	P-745	



MVR

Gauge	Bevel	Blade	Reusable
20g	Double	Straight	07520
23g	Double	Straight	07523



Multisided Scleral

Bevel	Thickness	Reusable
Double	0.38mm	0780



Paracentesis

Width	Bevel	Blade	Reusable
1.0mm	Up	Straight	07610
1.0mm	Up	Angled	07610BU



PREMIER EDGE® MICROSURGICAL KNIVES

Newly redesigned Premier Edge® microsurgical safety knives support compliance with OSHA 1910.1030(d)(2)(i) and 1910.1030(d)(2)(vii)(B), and are CMS compliant. The scalpels have a patented and ergonomically designed safety shield with ribbed thumb-grip. The natural sliding thumb motion allows one-handed operation and ambidextrous usability. The safety shield secures with a distinctive click.

Disposable	■ 6 per box	■ Single use	■ Sterile	
Safety	■ 6 per box	■ Single use	■ Sterile	■ One-handed retractable safety guard

PATENT NO. 8,875,405 PATENT NO. D639,432



Slit

Width	Bevel	Blade	Disposable
1.8mm	Up	Angled	PE3818
2.2mm	Up	Angled	PE3822
2.4mm	Up	Angled	PE3824
2.5mm	Up	Angled	PE3825
2.65mm	Up	Angled	PE3826
2.75mm	Up	Angled	PE3827
2.85mm	Up	Angled	PE3828
3.0mm	Up	Angled	PE3830
3.2mm	Up	Angled	PE3832
3.2mm	Up	Straight	PE3832S



Clear Cornea Slit w/ Depth Marking

Width	Bevel	Blade	Disposable	Safety
2.2mm	Double	Angled	PE3122	
2.5mm	Double	Angled	PE3125	
2.65mm	Double	Angled	PE3126	PE3126-SL
2.75mm	Double	Angled	PE3127	PE3127-SL
2.85mm	Double	Angled	PE3128	PE3128-SL
3.2mm	Double	Angled	PE3132	



Clear Cornea Parallel Slit

Width	Bevel	Blade	Disposable
2.65mm	Double	Angled	PE3226
2.75mm	Double	Angled	PE3227
2.85mm	Double	Angled	PE3228
3.0mm	Double	Angled	PE3230



Thin Slit

Width	Bevel	Blade	Disposable	Safety
2.0mm	Up	Angled	PE4820-TL	PE4820-SL
2.2mm	Up	Angled	PE4822-TL	PE4822-SL
2.4mm	Up	Angled	PE4824-TL	PE4824-SL
2.65mm	Up	Angled	PE4826-TL	PE4826-SL
2.75mm	Up	Angled	PE4827-TL	PE4827-SL
2.85mm	Up	Angled	PE4828-TL	PE4828-SL
3.0mm	Up	Angled	PE4830-TL	PE4830-SL
3.2mm	Up	Angled	PE4832-TL	PE4832-SL



Clear Cornea Trapezoid Slit

Width Mark	Bevel	Blade	Disposable
1.5mm - 1.7mm	Up	Angled	PE3303
1.8mm - 2.2mm	Up	Angled	PE3301
2.5mm - 3.0mm	Up	Angled	PE3302



Implant

Width	Bevel	Blade	Disposable
4.0mm	Up	Angled	PE3440



Round Tunnel

Width	Bevel	Blade	Disposable
2.2mm	Double	Angled	PE3622
3.0mm	Double	Angled	PE3630



Stab

Angle	Bevel	Blade	Disposable	Safety
15°	Double	Straight	PE3015	PE3015-SL
22.5°	Double	Straight	PE3022	
30°	Double	Straight	PE3030	
45°	Double	Straight	PE3045	



MVR

Gauge	Bevel	Blade	Disposable	Safety
19g	Double	Straight	PE3519	
20g	Double	Straight	PE3520	PE3520-SL
23g	Double	Straight	PE3523	
25g	Double	Straight	PE3525	
19g	Double	Angled	PE3519A	
20g	Double	Angled	PE3520A	
23g	Double	Angled	PE3523A	



Restricted Depth Stab

Angle	Depth	Disposable
15°	3.0mm	PE3015-3
15°	5.0mm	PE3015-5



Side Port



Width	Bevel	Blade	Disposable
1.0mm	Double	Angled	PE3010
1.2mm	Double	Angled	PE3012
1.4mm	Double	Angled	PE3014







Crescent Tunnel

Width	Bevel	Blade	Disposable	Safety
2.0mm	Up	Angled	PE3720A	PE3720A-SL
2.0mm	Up	Straight	PE3720	

Myringotomy Blades

Description	Ref #
 Spear Tip	PE7120
 Lance Tip	PE7121

Mini Blades

Description	Ref #
 Multi-Sided, sharp all around	PE5700
 Rounded Tip, sharp on one side	PE6400
 Curved Tip, sharp on one side	PE6700
 Rounded Tip, sharp all around	PE6900



NEEDLES, CANNULAS & CYSTOTOMES

■ 10 per box

■ Single use

■ Sterile

Anesthesia

	Description	Gauge	Length	Ref #
	Retrobulbar	25g (0.50)	1 1/2" (38mm)	4004
	Retrobulbar Atkinson Tip	25g (0.50)	1 1/2" (38mm)	4004B
	Retrobulbar Atkinson Tip	23g (0.60)	1 1/2" (38mm)	4004C
	Peribulbar Atkinson Tip	25g (0.50)	7/8" (22mm)	4004D
	Peribulbar Atkinson Tip	23g (0.60)	1 1/4" (32mm)	4004F
	Peribulbar Atkinson Tip	25g (0.50)	1 1/4" (32mm)	4004M
	Sub Tenon	19g (1.10)	1" (25mm)	4005
	Tri-Port Sub Tenon	19g (1.10)	1" (25mm)	4005TP

Irrigation/Air Injection & Viscoelastic

	Description	Gauge	Length	Ref #
	Anterior Chamber	25g (0.50)	7/8" (22mm) 4mm bend	4001AF
	Anterior Chamber	27g (0.40)	7/8" (22mm) 4mm bend	4002F
	Anterior Chamber	30g (0.30)	7/8" (22mm) 4mm bend	4001F
	Anterior Chamber	30g (0.30)	7/8" (22mm) 6mm bend	4001
	Anterior Chamber	25g (0.50)	7/8" (22mm) 6mm bend	4001A
	Anterior Chamber	27g (0.40)	7/8" (22mm) 6mm bend	4002
	Anterior Chamber/Viscoelastic Thin Wall	23g (0.60)	7/8" (22mm) 8mm bend	4001CTW
	Anterior Chamber/Viscoelastic	27g (0.40)	7/8" (22mm) 8mm bend	4002B
	Viscoelastic	27g (0.40)	7/8" (22mm) 9mm bend	4002L
	Anterior Chamber/Viscoelastic	23g (0.60)	7/8" (22mm) 10mm bend	4001B
	Anterior Chamber/Viscoelastic	27g (0.40)	7/8" (22mm) 11mm bend	4002A
	Curved	25g (0.50)	1 1/8" (29mm)	4003
	Anterior Chamber Washout	19g (1.10)	1" (25mm) 8mm bend	4036

Infusion Cannula

	Description	Ref #
	2.0mm length	4042-2S
	4.0mm length	4042-4S

NEEDLES, CANNULAS & CYSTOTOMES










■ 10 per box









■ Single use

■ Sterile

Capsulotomy

	Description	Gauge	Length	Ref #
	Straight Irrigating	25g (0.50)	7/8" (22mm)	4017A
	Straight Irrigating	27g (0.40)	5/8" (16mm)	4017IS
	Formed Irrigating	25g (0.50)	7/8" (22mm)	4017
	Formed Irrigating	27g (0.40)	7/8" (22mm)	4017B
	Formed Irrigating	25g (0.50)	5/8" (16mm)	4017G
	Formed Irrigating	27g (0.40)	5/8" (16mm)	4017IO
	Formed Irrigating	30g (0.30)	5/8" (16mm)	4017I
	Formed Irrigating	25g (0.50)	5/8" (16mm) Radius of curvature 1/2"	4017IE
	Angled Irrigating	25g (0.50)	5/8" (16mm) 12mm bend	4017M
	Angled Irrigating	27g (0.40)	5/8" (16mm) 12mm bend	4017N
	F-SL Curve Irrigating	25g (0.50)	5/8" (16mm) 19mm radius	4017Z
	Reverse Formed Irrigating	25g (0.50)	5/8" (16mm)	4017L
	Formed Capsulorhexis	25g (0.50)	5/8" (16mm)	4017S
	Formed Capsulorhexis	27g (0.40)	5/8" (16mm)	4017R

Hydrodissection

	Description	Gauge	Length	Ref #
	U-Shaped	25g (0.50)	7/8" (22mm)	4016AH
	Horizontal Flow, flattened tip	25g (0.50)	7/8" (22mm) 8mm bend	4036A
	Horizontal Flow, flattened tip	27g (0.40)	7/8" (22mm) 8mm bend	4036B
	Helsinki, curved with flattened tip (Horizontal Flow)	25g (0.50)	7/8" (22mm)	4036F
	Helsinki, curved with flattened tip (Horizontal Flow)	27g (0.40)	7/8" (22mm)	4036G
	Vertical Flow, flattened tip	25g (0.50)	7/8" (22mm) 4mm bend	4036C
	Vertical Flow, flattened tip	27g (0.40)	7/8" (22mm) 4mm bend	4036D
	Sauter, curved with flattened tip (Vertical Flow)	25g (0.50)	7/8" (22mm)	4036M
	Sauter, curved with flattened tip (Vertical Flow)	27g (0.40)	7/8" (22mm)	4036N
	Cortical Cleaving, round end	27g (0.40)	7/8" (22mm) 11mm bend	4036P
	Cortical Cleaving, flat end	27g (0.40)	7/8" (22mm) 11mm bend	4036Q
	Kellan, curved, 20mm radius, flat 3mm from end	25g (0.50)	7/8" (22mm)	4068
	Chang, 90° angled flattened tip	27g (0.40)	3/4" (19mm) 1mm bend	4036J









NEEDLES, CANNULAS & CYSTOTOMES

■ 10 per box





■ Single use

■ Sterile


Capsule Polishers

	Description	Gauge	Length	Ref #
	Kratz, 0.3mm top port	21g (0.80)	7/8" (22mm) 9mm bend	4027
	Kratz, 0.3mm top port	23g (0.60)	7/8" (22mm) 9mm bend	4027A
	Kratz, open end	21g (0.80)	7/8" (22mm) 6mm bend	4027B
	Olive-Tip Angled	25g (0.50)	1" (25mm) 10mm bend	4027E
	Terry	27g (0.40)	7/8" (22mm) 6mm bend	4028
	Capsule Polisher	23g (0.60)	7/8" (22mm) 7mm bend	4028B


Cortex Extractors

	Description	Gauge	Length	Ref #
	Simcoe	23g (0.60)	1/2" (13mm)	4008
	U-Shaped	21g (0.80)	7/8" (22mm)	4016
	U-Shaped	23g (0.60)	7/8" (22mm)	4016F
	U-Shaped	27g (0.40)	7/8" (22mm)	4016B
	Micro J-shaped	25g (0.50)	7/8" (22mm)	4016G
	Charleux	23g (0.60)	1/2" (13mm)	4026


Anterior Chamber Maintainer

	Description	Gauge	Length	Ref #
	Anterior Chamber Maintainer	20g (0.90)	8.5" silicone tubing with 3mm tip	4300P

Irrigating & Aspirating

	Description	Gauge	Ref #
	Simcoe Double I & A Cannula, 0.30mm aspiration port	23g/23g (0.60/0.60)	4007

Lens Loops

	Description	Gauge	Length	Ref #
	Irrigating Vectis, three 0.25mm ports	24g (0.55)	5.5mm	4029A













LASIK Cannulas

■ 10 per box

■ Single use

■ Sterile

	Description	Gauge	Length	Ref #	
	Lin, 9mm flattened angle, two 0.3mm ports	25g (0.50)	7/8" (22mm)	4036L	
	Buratto, 11mm angle, one front and two side ports	25g (0.50)	7/8" (22mm)	4036LB	
	Lindstrom, 45° angle, 5/8" radius, three 0.3mm ports	25g (0.50)	7/8" (22mm)	4036LR	
	Slade, 9mm formed extension with spatulated tip	26g (0.45)	1" (25mm)	4036LS	
	Open End, six 0.15mm side ports	25g (0.50)	7/8" (22mm)	4225	

SOFT CELL® PLUS and SOFT CELL® PVA Foam

NEW! SOFT CELL PLUS is a sterile PVA Foam. Fast wicking with dense form, fast absorbency with low particulates.

SOFT CELL PVA Foam is sterile, fast absorbency with low particulates.



Qty	Description	Ref #
150 per box	SOFT CELL Surgical Spears	0525
150 per box	SOFT CELL PLUS Surgical Spears	0530
20 per box	SOFT CELL PLUS Wipe + Wick	0571
20 per box	SOFT CELL PLUS LASIK Eye Drain (Chayet) 10mm internal diameter	0572
20 per box	SOFT CELL PLUS LASIK Eye Drain (Chayet/Katzen) 12.2mm internal diameter	0573
200 per box	SOFT CELL PLUS Eye Pledget	0574
20 per box	SOFT CELL PLUS LASIK Shield, 9mm	0575

VISCO SHIELD® Viscoelastic

FOR EXPORT ONLY



VISCO SHIELD is a dispersive HPMC viscoelastic composed of a sterile, high molecular weight, hydroxypropyl methylcellulose solution for use in anterior segment surgical procedures including cataract extraction and intraocular lens implantation.

- 2ml filled syringe with cannula
- HPMC (hydroxypropyl methylcellulose)
- Stores at room temperature

Qty	Description	Ref #
6 per box	20,000 cps	5122
6 per box	8,000 cps	51082

VISCO SHIELD® Topical

FOR EXPORT ONLY



VISCO SHIELD Topical is preservative-free and designed to keep the cornea hydrated. It adheres to the cornea and forms a thin, optically neutral clear film.

- 1.5ml filled syringe with cannula
- HPMC (hydroxypropyl methylcellulose)
- Stores at room temperature

Qty	Description	Ref #
6 per box	8,000 cps	5200

Surgical Accessories

Single use.



Qty	Description	Ref #
30 per box	Iris Retractors	4465
10 per box	OASIS Eye Cautery Fine tip. Battery operated. Low temperature. Range: 1100°F ± 50°F	3501-A
10 per box	Pre-inked Corneal Marking Pen	0442
50 per box	Disc Fluid Filter Filters out particulate and biological contaminate matter larger than 0.22 µ.	0601

Marketing Materials



Description	Ref #
Dry Eye & Surgery Patient Brochure	LIT-6540
Patient Brochure, Intacs English	LIT-INTACS-PB
Patient Brochure, Intacs Spanish	LIT-INTACS-PBSP
Product Sheet for MDs, Intacs	LIT-INTACSM-D-S
Product Sheet for ODs, Intacs	LIT-INTACSOD-S
Wall Chart, OASIS Iris Expander	LIT-WCIE1013



Donawa Lifescience Consulting Srl
Piazza Albania 10
00153 Rome, Italy



Glendora, CA, USA

(800) 528-9786 (USA Toll Free) | (909) 305-5400

customerservice@oasismedical.com | www.oasismedical.com



EN ISO 13485:2012 + AC : 2012
LIT-SRG 04/15