

Patient chair Malin 4



GENERAL

Malin 4 patient chair is popular in several treatment areas because it combines good ergonomics and a high functionality. The stable design with large wheels, robust mechanical design provides safety for both patient and staff.

Malin 4 is electric and easy to adjust in height with the hand control. Malin 4 is available in different models and with plenty of accessories to tailor make the chair. The chair has a manual brake bar on the back wheels, electric brake is available as an option.

MALIN 4 ROCK

120-00103-00



Malin 4 Rock is our appreciated base model that can be outfitted with accessories for different examination situations. Seat and back rest can be tilted to get the desired angle.

Lifting capacity	150 kg / 331lb
Weight	41kg / 90lb
Height adjustment	540-840mm / 21-33in
Wheel frame	530x580mm / 21x23in
Backrest angle	-8° to +20°
Seat angle	+4° to -10°
Battery operated	

MALIN 4 SWING

120-00103-01



Malin 4 Swing has a swivel seat that can be rotated $\pm 90^\circ$, to enable the patient to comfortably be seated from the side with the chair parked. The seat is conveniently turned to the correct position by the device table while the base of the chair remains in position.

Lifting capacity	150 kg / 331lb
Weight	41kg / 90lb
Height adjustment	540-840mm / 21-33in
Wheel frame	530x580mm / 21x23in
Battery operated	

MALIN 4 SLIM

120-0103-02



Malin 4 Slim has a narrower wheelbase in the front which makes it easier to manoeuvre the chair on places with limited space, e.g. close to device tables with only one lifting column.

Lifting capacity	150 kg / 331lb
Weight	41kg / 90lb
Height adjustment	540-840mm / 21-33in
Wheel frame	530x580mm / 21x23in
Backrest angle	-8° to +20°
Seat angle	+4° to -10°
Battery operated	

STANDARD SEAT COLOURS



140-00135-00 blue vinyl

Vinyl has an airy, ventilated slick texture.

Seat 450x470mm / 18x19in, Backrest 360x410mm / 14x16in



140-00247-02 black art leather

Artificial leather is more durable.

ACCESSORIES

There are a variety of accessories that all have been developed in close cooperation with medical professions, to optimize the chair for different treatments, to get an ergonomic and comfortable environment for both patients and staff.

140-00128-00 Extended lifting height, 100 mm / 4in extra

140-00126-02 Wireless hand control

140-00126-03 Wireless Hand control for electric brake

IR hand control comes with a strap and bracket/holder to be placed where its suitable, on the chair or by the working station.

140-00250-00 Electric brake

One button brake and release button on either the fixed or IR hand control

140-00315-00 Wheels with direction control

Makes the chair easier to move around when back wheels have a direction lock



150-00102-00 Armrest standard left

150-00102-01 Armrest standard right

150-00103-00 Armrest collapsible left

150-00103-01 Armrest collapsible right



150-00105-00 Armrest for fluorescein photography left

150-00105-01 Armrest for fluorescein photography right



140-00129-03 Foot support, Rubber (not for Malin slim)

The foot support raises with the chairs and can be adjusted in height.




150-00132-00 Driving handle 450mm

150-00132-02 Driving handle 600mm



190-00018-00 Head support



rini

Rini Ergoteknik AB
Fågelsångsvägen 3B
SE-186 42 Vallentuna
Sweden
Phone +46 8 594 77170
Email info@rini.se
www.rinicompany.com



rini
MedTech

Rini Medtech Inc
330 Summerview Ct San Ramon
CA94583
USA
Phone +1-925-905-1320
Email info@rinimedtech.com
www.rinimedtech.com

