



rini

Medical Ergonomics made in Sweden



www.rinicompany.com

Product catalogue 2024
issue 3



| Table Of Contents | Side |
|-----------------------------------|------|
| 1. Operating tables - RiEye Mk2S | 3 |
| 2. Surgeon chairs - Carl | 5 |
| 3. Patient chairs & Device tables | 7 |
| 4. Other MedTech-products | 9 |
| 5. General Terms and Conditions | 13 |

Welcome to Rini!

We are a Swedish company with innovative research and development and have since mid-1990 been a quality supplier of ergonomic products and solutions to medical professionals. We have our own sales and support organization in Sweden and work mainly with appointed distributors on our export markets.

Our business strategy is simple!

Rini offers ergonomic products and solutions uniquely tailored to different treatment areas. Giving doctors and associated personnel a maximum comfortable and productive work environment.

The close contact we have with many highly qualified users is exciting and a positive challenge to continue to develop our products further - always with customer satisfaction in focus! **Read more at www.rinicompany.com**

Rini is a quality certified supplier according to ISO 13485.



Contact us by phone +46 8 594 77170 or E-mail info@rini.se and we will tell you more!

1. Operating tables

Operating table RiEye Mk2S

Product code Description

- 100-00105-05 RiEye Mk2S Operating table R5 – Suitable for shorter treatments such as injections
- Vertical leg support make it easier for the patients to take a seat in the chair straight from the front
 - Quick and safe positioning of the patient with Rini's patented "HandsFree" head support
 - Headrest with a unique ergonomic design
 - Option to save 6 user positions
 - Steady takes patient weights up to 300Kg
 - Adaptive chassis means RiEye Mk2S stands firmly on uneven floors
 - The chairs are mobile and can also be moved sideways
 - Short leg cushion 400mm long



- 100-00105-06 RiEye Mk2S Operating table – R6 Suitable for longer treatments
- Slightly angled leg pad in the lowest position, patient sits down from the front or side
 - Quick and safe positioning of the patient with Rini's patented "HandsFree" head support
 - Headrest with a unique ergonomic design
 - Option to save 6 user positions
 - Steady takes patient weights up to 300Kg
 - Adaptive chassis means RiEye Mk2S stands firmly on uneven floors
 - The chairs are mobile and can also be moved sideways
 - Long leg cushion 500mm long



- 100-00105-07 RiEye Mk2S Operating table R7 – Suitable for longer treatments, including anesthesia operations
- Angled leg pad in lowest position patient sits down in the chair from side
 - Quick and safe positioning of the patient with Rini's patented "HandsFree" head support
 - Headrest with a unique ergonomic design
 - Option to save 6 user positions
 - Steady takes patient weights up to 300Kg
 - Adaptive chassis means RiEye Mk2S stands firmly on uneven floors
 - The chairs are mobile and can also be moved sideways
 - Extra long leg cushion 680mm



Colors selection for ILLUSION cushion kits, Phthalate-free PU, Martindale 500.000

| Product code | Description | R5 | R6 | R7 |
|--------------|-----------------------------------|-----|-----|-----|
| 140-00331-xx | ILLUSION cushion kit blue | -00 | -01 | -02 |
| 140-00333-xx | ILLUSION cushion kit black | -00 | -01 | -02 |
| 140-0033x-xx | ILLUSION cushion kit other colors | -00 | -01 | -02 |
| 140-00285-xx | Extra wide cushion kit | -00 | -01 | -02 |

Accessories for operating table RiEye Mk2S

Accessories marked with * must be factory mounted

Product code Description

Headrest and headrest accessories

| | |
|--------------|--|
| 140-00270-00 | Headrest X-Tray Doublelink, replaces Bowl (standard) headrest |
| 140-00270-01 | Pillow for X-Tray black (300x265mm WxL) |
| 140-00270-02 | Pillow for X-Tray grey |
| 140-00270-03 | Headrest S-Tray Doublelink, replaces Bowl (standard) headrest |
| 140-00270-04 | Pillow for S-Tray (210x180mm WxL) |
| 140-00270-14 | New Pillow for S-tray, U-shaped (220x220mm) |
| 140-00252-02 | Headrest S-Tray Nippon, replaces Bowl (standard) headrest |
| 140-00252-05 | Headrest FaceDown/FaceUp Doublelink, replaces Bowl (standard) headrest |
| 140-00270-06 | Pillow for FaceDown/FaceUp |
| 140-00270-07 | New Headrest OneGrip, replaces Doublelink Bowl(standard) headrest |
| 140-00270-08 | New Cushion Luigi for Onegrip headrest |
| 140-00270-09 | New S-tray for onegrip headrest |

| | |
|--------------|--|
| 140-00229-01 | Electrical headrest for Doublelink* |
| 140-03420-00 | Bowl headrest (standard) with Foundation for headrests to make headrest dockable |
| 140-00235-02 | Accessory rail for X-Tray (for wrist support) |
| 140-00242-01 | Wrist support Bowl(standard) / S-Tray for Doublelink |
| 140-00242-02 | Wrist support X-Tray |

Hand control holder

| | |
|--------------|--|
| 140-00238-01 | Holder for hand control unit flexible, bracket included |
| 140-00238-02 | Holder for hand control unit fixed, bracket included |
| 140-00238-03 | Holder for hand control unit fixed for IV pole, bracket included |

Foot control

| | |
|--------------|---|
| 140-00230-01 | Foot control with two joysticks movable tray |
| 140-00230-02 | Foot control with two joysticks fixed tray, cannot be combined with electrical headrest |

Personal comfort

| | |
|--------------|--|
| 140-00261-01 | Auto access to memory location |
| 140-00266-01 | Electrical brake and direction control * |
| 140-00352-00 | Reverse trendelenburg max 12° * |

Patient comfort

| | |
|--------------|---------------------------------------|
| 140-00350-00 | Width extension back rest right 410mm |
| 140-00350-01 | Width extension back rest left 410mm |
| 140-00232-00 | Length extension back rest 120mm |
| 140-00265-00 | Leg support circular ø 150mm |
| 140-00265-01 | Leg support half circular 200mm high |

Transport of the table

| | |
|--------------|---|
| 140-00259-00 | Driving handle backrest foldable R/L |
| 140-00256-00 | Driving handle leg foldable R/L R6 and R7 |
| 140-00256-01 | Driving handle leg R/L R5 |

General accessories

| | |
|--------------|--|
| 140-00236-02 | Drape support, formable flexible |
| 140-00236-03 | Drape support, formable with air tube connection |
| 140-00235-00 | Attachment rail for accessories on the seat R or L 250mm |
| 140-00235-01 | Attachment rail for accessories on the seat R or L 250mm if Side rail is mounted |

| | |
|--------------|--|
| 140-00237-02 | Battery lead (1pc is standard) |
| 140-00237-04 | Battery Lithium iron, low weight and higher capacity |

| | |
|--------------|---|
| 140-00351-00 | Bracket for oxygen tank |
| 140-00163-00 | Anaesthesia armrest, clamp and arm strap included |
| 140-00239-01 | IV pole, bracket included |
| 140-00251-00 | Leg cushion cover R6 R7 |
| 140-00251-01 | Leg cushion cover R5 |

Patient safety

| | |
|--------------|---|
| 140-00240-00 | Side rail collapsible right side |
| 140-00240-01 | Side rail collapsible left side |
| 140-00320-00 | Armrest for Side rail collapsible * |
| 140-00240-02 | Side rail right or left side 600mm, 140-00235-00 included |
| 140-00263-00 | Belt waist |
| 140-00263-01 | Belt arm 700mm |
| 140-00263-02 | Belt leg 2100mm |

2. Surgeons chairs

Surgeon chairs Carl MK2, a new generation of electrically controlled chairs

Product code Description

- 120-00133-00 Carl Mk2 R6 - Rilis arm support
- Uniquely designed frame with three wheels for high availability
 - Electric controlled central brake
 - Seat height from 490 to 800mm
 - Seat 420x460mm back support 300x350mm
 - Battery-operated, charger is connected to chair
 - Back and seat can be adjusted in height and angle
 - The chair is equipped with Rilis safety arm support



Product code Description

- 120-00133-01 Carl Mk2 R7 - One hand grip support
- Uniquely designed frame with three wheels for high availability
 - Electric controlled central brake
 - Seat height from 490 to 800mm
 - Seat 420x460mm back support 300x350mm
 - Battery-operated, charger is connected to chair
 - Back and seat can be adjusted in height and angle
 - The chair is equipped with one hand grip support



Product code Description

- 120-00133-02 Carl Mk2 R8 - FlexiDoc
- Uniquely designed frame with three wheels for high availability
 - Electric controlled central brake
 - Seat height from 490 to 800mm
 - Seat 420x460mm back support 300x350mm
 - Battery-operated, charger is connected to chair
 - Back and seat can be adjusted in height and angle
 - With FlexiDoc armrests for the most demanding procedures



Product code Description

- 120-00133-03 Carl Mk2 R5 - Easy arm support
- Uniquely designed frame with three wheels for high availability
 - Electric controlled central brake
 - Seat height from 490 to 800mm
 - Seat 420x460mm back support 300x350mm (LxW)
 - Battery-operated, charger is connected to chair
 - Back and seat can be adjusted in height and angle
 - The chair is equipped with Easy arm support



Upholstery and color (state when ordered)

| Product code | Description | Blue | Black |
|--------------|--------------------------|-------------|-------|
| 140-00118-xx | Cushion kit vinyl | -00 | -02 |
| 140-00247-xx | Cushion kit art leather | -00 | -02 |
| 140-00306-xx | Cushion kit soft padding | Not offered | -02 |
| 140-00321-xx | Saddle seat Amazone | Not offered | -02 |
| 140-00322-xx | Saddle seat Jumper | Not offered | -02 |

Vinyl has an airy, ventilated slick texture, good for long sessions

Artificial leather has a denser structure and is easier to clean

Options for One hand grip armrest R7

Product code Description

- 140-00323-00 Arm cushion soft and convex shape
 140-00323-01 Arm cushion soft and rectangular shape 400x110mm (LxW)
 140-00324-00 Extended height adjustment 350mm (225mm)

Option for Carl Mk2

Product code Description

- 140-00327-00 Pedal platform 600x680mm (WxL) can be retrofitted

Surgeon chairs Carl 4 - electrically controlled chairs

| Product code | Description |
|--------------|---|
| 120-00102-04 | Carl 4 Foot R5 - Easy arm support |
| 120-00102-00 | Carl 4 Foot R6 - Rilis arm support |
| 120-00102-02 | Carl 4 Foot R7 - One hand grip support |
| 120-00102-05 | Carl 4 Heel R5 - Easy arm support |
| 120-00102-01 | Carl 4 Heel R6 - Rilis arm support |
| 120-00102-03 | Carl 4 Heel R7 - One hand grip support |
| | <ul style="list-style-type: none"> • Low friction, large 100mm castors • Seat height from 490 to 800mm • Seat 420x460mm back support 300x350mm (LxW) • Back and seat can be adjusted in height and angle • Battery-operated, charger is connected to chair • Mechanical brake, on the back wheels |



Foot R7



Heel R6



Heel R5

Surgeon chairs Carl - manually controlled chairs

| Product code | Description |
|--------------|--|
| 120-00100-00 | Carl Spring R6 - Rilis arm support |
| 120-00100-01 | Carl Spring R7 - One hand grip support |
| | <ul style="list-style-type: none"> • Foot controlled height adjustment with ring • Seat height from 540 to 680mm • Seat 420x460mm back support 300x350mm • Seat and back can be adjusted in height and angle • Mechanical brake, on the back wheels |



R6



R7

| | |
|--------------|---|
| 120-00101-00 | Carl Swing R6 - Rilis arm support |
| 120-00101-01 | Carl Swing R7 - One hand grip support |
| | <ul style="list-style-type: none"> • Foot controlled height adjustment with ring • Seat height from 500 to 700mm • Seat 420x460mm back support 300x350mm • An easily movable and rotatable chair • Back and seat can be adjusted in height and angle |



R6



R7

Upholstery and color (state when ordered)

| Product code | Description | Blue | Black | |
|--------------|--------------------------|-------------|-------|----------------------------|
| 140-00118-xx | Cushion kit vinyl | -00 | -02 | |
| 140-00247-xx | Cushion kit art leather | -00 | -02 | |
| 140-00306-xx | Cushion kit soft padding | Not offered | -02 | (For Carl 4 Foot och Heel) |
| 140-00321-xx | Saddle seat Amazone | Not offered | -02 | (For Carl 4 Foot och Heel) |
| 140-00322-xx | Saddle seat Jumper | Not offered | -02 | (For Carl 4 Foot och Heel) |

Vinyl has an airy, ventilated slick texture, good for long sessions
Artificial leather has a denser structure and are easier to clean

140-00129-01 **New** Footrest Carl 4 Foot

Options for One hand grip armrest R7

| Product code | Description |
|--------------|--|
| 140-00323-00 | Arm cushion soft and convex shape |
| 140-00323-01 | Arm cushion soft and rectangular shape 400x110mm (LxW) |
| 140-00324-00 | Extended height adjustment 350mm (225mm) |

3. Patient chairs & Device tables

Patient chairs Malin 4 - electrically controlled

Product code Description

- 120-00103-00 Malin 4 Rock - hand controlled
- Low friction large 100mm castors with double ball bearing wheels
 - Lifting capacity from 540 to 840mm
 - Seat 450x470mm back support 360x410mm (LxW)
 - Battery powered
 - Stable manage patient weights up to 150Kg



- 120-00103-01 Malin 4 Swing - hand controlled
- Seat can be rotated $\pm 90^\circ$
 - Low friction large 100mm castors with double ball bearing wheels
 - Lifting capacity from 540 to 840mm
 - Seat 450x470mm back support 360x410mm
 - Battery powered
 - Stable manage patient weights up to 150Kg



- 120-00103-02 Malin 4 Slim - hand controlled
- Wheelbase for narrow floor space under the table
 - Low friction large 100mm castors with double ball bearing wheels
 - Lifting capacity from 540 to 840mm
 - Seat 450x470mm back support 360x410mm
 - Battery powered
 - Stable manage patient weights up to 150Kg



Options:

Product code Description

- 140-00128-01 Extended lifting height, 100 mm extra height adjustment
- 140-00126-02 Wireless hand control infrared
- 140-00126-03 Wireless hand control infrared, if electrical brake
- 140-00250-01 Electrical brake by the standard hand control
- 140-00315-00 Wheels with direction control

Accessories:

Product code Description

- 150-00102-00 Armrest standard left
- 150-00102-01 Armrest standard right
- 150-00103-00 Armrest collapsible left
- 150-00103-01 Armrest collapsible right
- 150-00105-00 Armrest for fluorescein photography left 350mm
- 150-00105-01 Armrest for fluorescein photography right 350mm

- 140-00129-03 Foot support plate of rubber, not for Malin Slim
- 140-00132-00 Driving handle 450mm
- 140-00132-02 Driving handle 600mm
- 140-00260-00 Belt for patient
- 190-00018-00 Head support



Colors of cushion kit:

Product code Description

- 140-00135-00 Cushion kit blue
- 140-00135-02 Cushion kit black

Device table Linnea
Product code Description

- 130-00103-00 Linnea 850 - examination table
- Width between columns 850mm
 - Table top of beech laminate, size 1000x420x22mm
 - Electrically adjustable height from 690 to 1040mm
 - Manage weight up to 80Kg
 - Connected to 230 V
 - Equipped with safety brake for the patient

- 130-00109-00 Linnea Piccolo - examination table,
- Width between columns 720mm
 - Table top of light grey laminate, size 870x370/420x22mm
 - Electrically adjustable height from 580 to 1240mm
 - Manage weight up to 50Kg
 - Connected to 230 V
 - Equipped with low wheel base and two sockets 230V outlets as standard


Examination table Linnea 850 table top variants:
Product code Description


- 140-00142-00 Table top of beech laminate 1000x350x22mm
- 140-00143-00 Table top of beech laminate 1000x500x22mm
- 140-00144-00 Table top of beech laminate 1180x350x22mm, asymmetric
- 140-00145-00 Table top of grey laminate with birch edge 1000x450x38mm

Examination table Linnea Piccolo table top variants:
Product code Description

- 140-00150-00 Table top of light grey laminate, size 870x420x22mm
- 140-00145-03 Table top of grey laminate with birch edge 1000x420x38mm cut out for H-S chin support

Options:
Product code Description

- 140-00139-00 Low wheel base, table lowered with 75mm
- 140-00121-00 Control box mounted to column for more space beneath table top
- 140-00248-00 Factory mounted power outlet under table top


Accessories:
Product code Description

- 140-00272-00 Chin support bracket H-S 18mm
- 140-00272-01 Chin support bracket H-S 60mm

4. Other MedTech-products

Special chairs for surgeons and assistants

Product code **Description**

- 120-00115-01 ComfortMove - foot controlled
 120-00115-00 ComfortMove - hand controlled
- Ideal for long seatings, popular among anaesthetists
 - Thick and comfortable padding
 - Seat 430x400mm back support 300x290mm (LxW)
 - Black artleather as standard, other as option

Standard castors

Castors 50mm low friction

| Gas spring | Low | Medium | High |
|---------------------|-----|--------|------|
| Max seat height, mm | 570 | 700 | 820 |
| Min seat height, mm | 430 | 500 | 580 |



- 120-00105-01 RiFlex - foot controlled
 120-00105-00 RiFlex - hand controlled
- Adjustable lower back support
 - Popular Assistant chair for eye or ear surgery
 - Seat 360x390mm back support 230x370mm

Standard castors

Castors 50mm low friction

| Gas spring | Low | Medium | High |
|---------------------|-----|--------|------|
| Max seat height, mm | 570 | 690 | 870 |
| Min seat height, mm | 430 | 500 | 580 |



| Description | Blue | Black | Other color |
|--------------------------------|------|-------|-------------|
| 140-00245-xx Cushion kit vinyl | -00 | -02 | - |

- 120-00104-01 Dynamic - foot controlled
 120-00104-00 Dynamic - hand controlled
- Chair with unique dynamic lower back support
 - Popular operator chair within dental care
 - Seat 400x440mm back support 160x140mm
 - Black stamskin as standard, other as option

Standard castors

Castors 50mm low friction, brake unloaded

| Gas spring | Low | Medium | High |
|---------------------|-----|--------|------|
| Max seat height, mm | 600 | 720 | 900 |
| Min seat height, mm | 460 | 530 | 620 |



- 120-00111-01 RiRest - foot controlled
 120-00111-00 RiRest - hand controlled
 120-00111-03 RiRest - backrest foot controlled
 120-00111-02 RiRest - backrest hand controlled
- Soft seat small and movable
 - Popular within operating theaters
 - Seat Ø330mm back support 230x370mm

Standard castors

Castors 50mm high friction

| Gas spring | Low | Medium | High |
|---------------------|-----|--------|------|
| Max seat height, mm | 600 | 720 | 900 |
| Min seat height, mm | 460 | 530 | 620 |



120-00111-01

120-00111-02

| Description | Blue | Black | Other color |
|--|------|-------|-------------|
| 140-00246-0x Cushion kit vinyl | -00 | -02 | - |
| 140-00246-1x Cushion kit vinyl with backrest | -10 | -12 | - |

Choice of castors**Product code Description**

140-00122-00 Castors 50mm low friction - set

140-00122-01 Castors 50mm high friction - set

Choice of gas spring**Product code Description**

140-00124-00 Gas spring low hand controlled chair

140-00124-01 Gas spring medium hand controlled chair

140-00124-02 Gas spring high hand controlled chair

140-00124-03 Gas spring low foot controlled chair

140-00124-04 Gas spring medium foot controlled chair

140-00124-05 Gas spring high foot controlled chair

Accessories**Product code Description**

140-00183-00 Foot support foot controlled chairs

140-00182-00 Foot support hand controlled chairs

Saddle chairs for surgeons and assistants

Product code Description

| | |
|--------------|--|
| 120-00106-01 | Amazone - saddle chair foot controlled |
| 120-00106-00 | Amazone - saddle chair hand controlled |
| 120-00106-02 | Amazone - back support hand controlled |
| 120-00106-03 | Amazone - back support foot controlled |
| 120-00106-04 | AmazoneBalance - hand controlled |
| 120-00106-05 | AmazoneBalance - foot controlled |

- those who finds other saddle chairs too wide
- Thick and comfortable padding and tiltable seat
- Seat 340x350mm back support 210x330mm (LxW)
- Black artleather as standard, other as option

Standard castors

Castors 65mm break unloaded

| Gas spring | Low | Medium | High |
|---------------------|-----|--------|------|
| Max seat height, mm | 570 | 690 | 820 |
| Min seat height, mm | 430 | 500 | 580 |



120-00106-03



120-00106-05

| | |
|--------------|---------------------------------------|
| 120-00110-01 | Jumper - saddle chair foot controlled |
| 120-00110-00 | Jumper - saddle chair hand controlled |
| 120-00110-02 | Jumper - back support hand controlled |
| 120-00110-03 | Jumper - back support foot controlled |
| 120-00110-04 | JumperBalance - hand controlled |
| 120-00110-05 | JumperBalance - foot controlled |

- A chair for longer sittings
- Excellent for those who wants a lower seat position
- Seat 360x430mm back support 210x330mm
- Black artleather as standard, other as option

Standard castors

Castors 65mm break unloaded

| Gas spring | Low | Medium | High |
|---------------------|-----|--------|------|
| Max seat height, mm | 570 | 690 | 820 |
| Min seat height, mm | 430 | 500 | 580 |



| | |
|--------------|---|
| 120-00123-01 | Perfect lite - saddle chair foot controlled |
| 120-00123-00 | Perfect lite - saddle chair hand controlled |

- Easy to maintain clean due to a seamless upholstery
- A small chair, excellent for those who want a high seat position
- Seat 300x390mm
- Black artleather as standard, other as option

Standard castors

Castors 50mm low friction

| Gas spring | Low | Medium | High |
|---------------------|-----|--------|------|
| Max seat height, mm | 570 | 690 | 820 |
| Min seat height, mm | 430 | 500 | 580 |



Choice of castors

Product code Description

| | |
|--------------|----------------------------------|
| 140-00122-00 | Castors 50mm low friction - set |
| 140-00122-01 | Castors 50mm high friction - set |

Choice of gas spring

Product code Description

| | |
|--------------|---|
| 140-00124-00 | Gas spring low hand controlled chair |
| 140-00124-01 | Gas spring medium hand controlled chair |
| 140-00124-02 | Gas spring high hand controlled chair |
| 140-00124-03 | Gas spring low foot controlled chair |
| 140-00124-04 | Gas spring medium foot controlled chair |
| 140-00124-05 | Gas spring high foot controlled chair |

Accessories

Product code Description

| | |
|--------------|-------------------------------------|
| 140-00183-00 | Foot support foot controlled chairs |
| 140-00182-00 | Foot support hand controlled chairs |

Assistant table Anna - for operating theatres

Product code Description

- 130-00105-00 Anna Tripod - assistant table with three wheels
- Height adjustment from 700 to 1100mm
 - Size of table top 700x500mm
 - Manage weight up to 20Kg
 - Individually braking castors, low friction
 - Battery powered, charger connect to table
 - Sound absorbing table top



- 130-00106-00 Anna Cross - assistant table with four wheels
- Height adjustment from 850 to 1350mm
 - Size of table top 1140x600 mm
 - Very stable, can manage up to 40Kg
 - Individually braking castors, low friction
 - Battery powered, charger connect to table
 - Sound absorbing table top



- 130-00106-01 Anna Cross - assistant table with fall edges and four wheels
- Fall edges so the Medical-guidewire for Endoscopic procedures remains on the table
 - Fall edges designed to prevent bacterial collections
 - Size of table top 1000x600 mm
 - Very stable, can manage up to 40Kg
 - Individually braking castors, low friction
 - Battery powered, charger connect to table
 - Sound absorbing table top



Assistant table variants:

Product code Description

- 140-00149-00 Height adjustment 850-1350mm table 130-00105-00
- 140-00151-00 Table top size 600x400mm table 130-00105-00
- 140-00152-00 Table top size 1300x500mm table 130-00106-00
- 140-00152-01 Table top size 1200x600mm table 130-00106-00

Armrests for operating theatre and examination

- Product code** **Description**
- 150-00100-00 Rilis armrest R6
- Unique safety armrest
 - Adjusted with great precision
 - Brackets fitting all of Rini's operator and assistants chairs



- 150-00113-00 One grip armrest R7
- Unique adjustment can be adjusted draped
 - Adjusted with great precision
 - Brackets fitting all of Rini's operator and assistants chairs



Armrest brackets

- Product code** **Description**
- 000-03067-00 Bracket right/left elliptical
- 000-01130-00 Bracket right/left rectangular 30x15mm
- 140-00303-00 Bracket right/left ComfortMove special chair
- 140-00303-01 Bracket R/L elliptical with rail
- 000-06302-01 Clove right side
- 000-06301-01 Clove left side

Undercarriages

- Product code** **Description**
- 000-03577-01 Under carriage saddle chairs elliptical
- Jumper, Perfect lite
- 140-00200-00 Under carriage RiFlex ,Rirest
- 000-01325-01 Under carriage saddle chairs 30x15mm
- Amazone without backrest



Armrest for operation and health care

Arm rests

- Product code** **Description**
- 000-03379-00 Armrest plate in molded seamless design
- 000-01327-00 Armrest plate arched
- 000-03395-00 Armrest plate 270x130mm
- 000-00431-00 Armrest plate 345x170mm

Arm

- Product code** **Description**
- 140-00127-00 Arm 75° 400mm, ball joint for curved armrest
- 140-00180-00 Arm 105° 400mm long



Bracket chair

- Product code** **Description**
- 000-03374-00 Bracket right elliptical
- 000-03375-00 Bracket left elliptical
- 000-03372-00 Bracket right rectangular 30x15mm
- 000-03373-00 Bracket left rectangular 30x15mm



Ladder Lars

Product code **Description**

- 160-00102-00 Step stool Lars with handle
- Height 300mm at highest step
 - Very stable, can manage up to 120Kg
 - Provided with anti-slip protection

- 160-00102-02 Step stool Lars 3 steps with handle
- Height 450mm at highest step
 - Very stable, can manage up to 120Kg
 - Provided with anti-slip protection



Product code **Description**

- 130-00108-00 Neptun bath trolley
- Height adjustment 700 till 1100mm
 - Individually braked wheels
 - Battery operated
 - Can manage up to 20Kg



5. General Terms and Conditions

1. Orders, Prices etc.

All prices are exclusive of VAT and other taxes or official duties. A quotation shall be binding for Rini for thirty (30) days from the date of the offer, unless otherwise stated in writing.

Any information on weight, measurement, capacity etc. in catalogues, prospects, folders are estimates only. All technical data is provided subject to design modifications.

A change of customs, taxes or other official duties, which results in increased cost to Rini, entitles Rini to increase the offered or agreed price.

2. Payment

Orders are to be pre-paid before shipment if not otherwise agreed.

For customers with an established credit line, payment is due strictly thirty (30) days from date of invoice. Any payments not made when due shall accrue interest from the due date until the date payment is made at a per annum rate equal to 18 percentage and will result in that future orders are be pre-paid.

3. Transportation and Packing

RINI handle packing and delivery of products according to "INCO TERMS Ex Works" of not otherwise agreed.

If other INCO TERMS are agreed for delivery, transportation, insurance and other handling fees will be charged in accordance with Rini's rates prevailing at the time of delivery. If the customer request express delivery, e.g. express cargo, courier, airfreight, the customer is obliged to compensate Rini for the actual cost of transportation in addition to the normal handling charge.

4. Risk of loss of Product

The risk of the product passes to the customer when Rini has delivered the product according to the agreed INCO TERM. Hence, Rini is only responsible for damage that may occur to the product up to the point in time when the risk passes to the customer.

5. Delivery time

The indicated delivery time starts when all terms in the purchase contract has been fulfilled. In case of pre-paid orders the delivery time starts when Rini has received payment. The product is considered to be delivered when risk of loss and damage passes to the customer in accordance with Section 4.

6. Inspection of Product

When a delivery is received the customer shall inspect the transportation packages and products for any transportation damages. If a transport damage has occurred, before the risk of the product has passed to the customer, the customer shall immediately report this to the shipping company and send a copy to Rini with full details of the damages including digital photo.

7. Warranty

RINI warrants that product sold at the time of the passing of risk as stipulated in Section 4 will be free from faults and shortcomings in construction, material and manufacture for 12 months. If the product is of the type requiring alignment or service, Rini is only responsible for faults remaining after such alignment or service has been carried out.

The customer shall immediately after the discovery of the fault or shortcoming make a claim to RINI in the form of a detailed description of the fault or shortcoming. However, such claim must be presented to Rini no later than 7 days after the day the customer discovered or ought to have discovered the fault or shortcoming.

After the receipt of a claim Rini will, at Rini's option, accept the return of the faulty product and reimburse the customer an amount equal to the purchase price, or repair or replace the faulty product.

Rini's responsibility for faults or shortcomings of delivered products is limited to the above. The customer is not entitled to any other remedies, but the above, in relation to faults or shortcomings.

8. Indirect Damages

The customer is not entitled to, and Rini is not in any event liable to pay, compensation for damages which the delivered product has caused to other property or to persons or any other consequential damages, including but not limited to loss of profit, loss of production or any other indirect damages.

In the event Rini would be held liable for compensation to third party for fault or safety shortcomings under the stipulations of, as it is called, product liability, the customer shall fully indemnify Rini for all damages finally awarded or settlement agreed, including Rini's costs in relation thereto, unless the damages caused by the fault or security shortcoming were due to gross negligence of Rini. Rini will give the customer notice of such claim and make it possible for the customer, if customer so wish, to influence the formulation of Rini's statement of defence or the terms to a possible settlement with the person presenting the claim, provided the customer issues a security, acceptable to Rini, for any costs or damages which can be imposed upon Rini.

9. Retention of Title

Products delivered shall remain the property of Rini until payment in full has been made.

10. Intellectual Property Rights

Customer is obliged to immediately notify Rini of any claim from third party alleging that products delivered under this Agreement infringes intellectual property rights (patent, registered design, trademark or copyright) of a third party. Provided customer present such claim immediately and provided customer grants Rini sole authority to defend or settle the case, Rini undertakes, at RINI expense, to defend the case, or come to a settlement, with the third party whose intellectual property rights has been infringed by the product delivered under this Agreement.

If it is held that a product delivered under this Agreement has infringed third party intellectual property rights Rini shall, at his own expense, and at Rini's option, obtain for the customer the right to continue using the product, replace or modify it so it becomes non-infringing or, in exchange for the return of the product, reimburse the price of the product. Rini is always entitled to deduct and the customer is obliged to pay an amount equal to the use of the product from the time of delivery to the claim was presented. If Rini has modified or developed a new product upon customer instructions or specifications, customer shall fully reimburse Rini for any and all costs and damages that might occur as a result of the modified or new product infringing third party intellectual property rights.

11. Force Majeure

If either party is prevented from performing under this Agreement due to force majeure the other party may not establish any sanctions for this reason. However, if the hindrance regards the customer's possibility to carry out payment Rini is entitled to withhold delivery of the product until payment has been effected. If the hindrance to effect payment last for more than four (4) weeks Rini is entitled to cancel the purchase.

The term "force majeure" means in this Agreement that the performance under the agreement is prevented, aggravated or delayed due to war, rebellion or riot, mobilization or unexpected call for military service of similar dimensions, requisition, seizure, exchange restrictions, export- or import restrictions, general scarcity of merchandise, shortage of transportation, labour conflict, restrictions on fuel and power, fire, faulty or delayed deliveries from subcontractors due to above causes, as well as any other circumstance beyond the control of the party.

12. Governing Law and Disputes

This agreement and any disputes hereunder shall be governed by Swedish law.

Any legal action relating to this Agreement and all legal relationships associated with it shall be decided by District Court of Stockholm as exclusive jurisdiction of first instance. However, RINI is entitled to turn to the district court or the Enforcement Service in the county or city where the customer has his legal domicile as regards claim for payment due.

The fact that a dispute has occurred does not exempt the customer to make payment due or to fulfil any other obligations under the agreement.

13. Miscellaneous

Either party's failure to enforce breach of this Agreement or to exercise its rights due to such breach shall not be deemed a waiver of such rights to enforce any future breach of agreement.

What is stated in this document as regards all matters contained is the full and exhaustive description of RINI's responsibilities and the rights of the customer. As regards such matters dealt with the customer is not entitled to any other remedy or compensation.

Stockholm May 2009