

# 122/TN Universal Trial Frame (NEO)

Universal trial frame with a renewed design, made of satin-finished metal and grey plastic. Even lighter and more functional, easy to use and compatible with 38 mm diameter lenses. The components in contact with the patient are soft and anatomical to ensure maximum comfort.



## FEATURES

Weight 52gr.

PD adjustment 50-80 mm

Axis scale division: 5°

Lens rotation 360°

Corneal vertex measurement range 0-14 mm

Total 5 lens holders: 3 front, 2 behind

Temple adjustable in extension / tilt

Nose rest adjustable in height / tilt with automatic position locking system

Supplied with 4 nose pads of different widths