


EASYOPHT 
true innovative solutions

the chain of excellence

C A T A L O G U E 2 5



Development, production and sale of medical devices, disposables and softwares worldwide.

EasyOpht is the division of **FE-Group Srl** specialized in the production of innovative devices.

The first product developed and patented is EasyTon Plus, a disposable prism of high quality for applanation tonometry.

Our business has been growing since 2006, year of foundation, along with our portfolio of instruments and solutions for Ophthalmologists and Opticians.

In 2016, FE-Group Srl was founded at first as a result of the union of Frastema Ophthalmics and EasyOpht, following the concept that union makes strength. Maintaining the personality of our different brands, we led the collaboration of Frastema Ophthalmics, EasyOpht, Easylac and FE-Service under a single 'hat' signed FE-Group Srl, aiming to offer different services from a single source.


Innovation, Assistance and Precision are our mantra and represent the daily goal of our staff. Thanks to our commitment and dedication, we maintain the Group Certifications ISO 13485, guarantee of quality.













| | |
|--|----|
| Easyton Plus | 6 |
| Fluorescein and Schirmer | 8 |
| EasyChart..... | 10 |
| MSD: | |
| <i>B trial lens plastic rim</i> | 14 |
| <i>Meniscus trial lens plastic rim</i> | 18 |
| <i>Single lens plastic rim</i> | 22 |
| <i>B trial lens metal rim</i> | 24 |
| <i>Accessories</i> | 26 |
| <i>Aspheric Lenses</i> | 30 |



COD. 100401 - EASYTON PLUS BOX 100 PCS

Disposable prism for applanation tonometry

-  **HIGH DEFINITION MEASURE**
-  **PRACTICAL**
-  **DISPOSABLE**
-  **SAFE**
-  **EASY**
-  **STERILE**



PRACTICAL

Everything you need for a complete sterile tonometry in one pack.

EASY

The tonometry prism fits easily on all types of Goldmann tonometers.

STERILE

Prevents risks of infections.

DISPOSABLE

No need to clean the tonometer prism. Simply throw it away. Your patients will be impressed with your attention to hygiene. You will be impressed with its speed of use.

SAFE

Sterilised with ethylene oxide, the tonometer prism is safe to use with no risk of allergic reaction or infection.


Each individual EASYTON blister contains a sterile tonometer prism tip. This is all that you need to perform an accurate and sterile tonometry procedure. For your economical convenience, EASYTON comes in an attractive dispenser holding 100 pieces.



1 Remove the tonometer prism from the sterile blister careful not to contaminate the tip.



2 Insert the tip into the tonometer with the line aligned to the zero position.



3 Perform the intraocular pressure measurement.



4 Dispose the used tonometer tip.

CE 0476

Patented

MADE IN ITALY 

EasyTest



COD. 700100 FLUORESCEIN | COD. 700101 SCHIRMER

Ophthalmic diagnostic strips

- **STERILE:**
EVERY PACK CONTAINS 100 STERILE SINGLE STRIPS
- **DISPOSABLE:**
DESIGNED TO BE USED ONLY ONCE FOR MAXIMUM SAFETY
- **PRACTICAL**
- **EASY TO USE**



DESCRIPTION

Fluorescein Sodium- 100 Ophthalmics Strips.

Each strip is impregnated with 1,0 mg of Fluorescein Sodium.

Strips for the evaluation of integrity of ocular surface and for contact lens fitting.

Schirmer Test - 100 Ophthalmics Strips.

Sterile ophthalmic strips for Dry Eye test.

PACKAGING

Each box contains 100 sterile strips.


COD. 700101



INDICATIONS:

- Used for measurement of tear production
- Graduated scale for an easy reading

COD. 700100



INDICATIONS:

- To examine the integrity of the corneal surface and for contact lens fitting

EasyChart




COD. 300152 - EASYCHART MEDICAL LIGHT

- High Resolution LCD 24" Monitor (1920×1080 24 bit)
- Visual angle: 178°
- Maximum brightness: 250 (cd/m²)
- Usable area: 520×290 mm
- IR Remote Controller
- Test Distance: from 1.5 to 6.5 mt
- Support VESA
- Logarithmic and Decimal System
- Eye Chart Tests with random letters
- Videos and demo images
- Software: Linux Operating System - MOS Software has the license to use it with the device

Optional accessories:
Pedestal cod. 300503

 **EasyChart**
medical
LIGHT



COD. 300113 - EASYCHART EVO

- High Resolution LCD 24" Monitor (1920×1080 24 bit)
- Visual angle: 178°
- Maximum brightness: 250 (cd/m²)
- Usable area: 520×290 mm
- IR Remote Controller
- Test Distance: from 1.5 to 6.5 mt
- Support VESA
- Logarithmic and Decimal Systems
- Custom sequences are possible
- Videos and demo images
- Test for Driving License renewal with exportable or printable report
- Eye Chart Tests with random letters and numbers
- Software: Linux Operating System - MOS Software has the license to use it with the device











Optional accessories:**Glare test kit cod. 300114****Pedestal cod. 300503**











**EasyChart
Evo**
**COD. 300150 / 300151- EASYCHART POLAR PLUS (LINEAR/CIRCULAR)**

- High Resolution LCD 24" Monitor (1920×1080 24 bit)
- Visual angle: 178°
- Maximum brightness: 250 (cd/m²)
- Usable area: 520×290 mm
- IR Remote Controller
- Test Distance: from 1.5 to 6.5 mt
- Support VESA
- Linear and circular polarization
- Logarithmic and Decimal Systems
- Custom sequences are possible
- Videos and demo images
- Haase Sequence (MKH)
- Eye Chart Tests with random letters and numbers
- Software: Linux Operating System - MOS Software has the license to use it with the device











Optional accessories:**Pedestal cod. 300503**

**EasyChart
Polar Plus**









| | | | | |
|---|---|---|---|---|
| Visual Acuity | ✓ | ✓ | ✓ |  |
| Single-line Visual Acuity | ✓ | ✓ | ✓ |  |
| Single-letter Visual Acuity | ✓ | ✓ | ✓ |  |
| ETDRS | ✓ | ✓ | ✓ |  |
| Kids Tests (images) | ✓ | ✓ | ✓ |  |
| Contrast sensitivity test Visual Acuity 0,125 | ✓ | ✓ | ✓ |  |
| Horizontal Chart | ✓ | ✓ | ✓ |  |
| Visual Acuity single Landolt C | ✓ | ✓ | ✓ |  |
| Visual Acuity Landolt C | ✓ | ✓ | ✓ |  |
| Single-chart Kids Test | ✓ | ✓ | ✓ |  |

| | | | | |
|--|---|---|---|---|
| Contrast Sensitivity Test Visual Acuity 0,3 | ✓ | ✓ | ✓ |  |
| Vertical Chart | ✓ | ✓ | ✓ |  |
| Visual Acuity Single Snellen E | ✓ | ✓ | ✓ |  |
| Visual Acuity Snellen E | ✓ | ✓ | ✓ |  |
| Kids Tests (shapes) | ✓ | ✓ | ✓ |  |
| Astigmatism Test | ✓ | ✓ | ✓ |  |
| Schober Test | ✓ | ✓ | ✓ |  |
| Bichromatic Test (letters, rings, E Snellen) | ✓ | ✓ | ✓ |  |
| Worth Lights (shapes, letters, numbers, images, Casale Test) | ✓ | ✓ | ✓ |  |
| Punctiform Chart (dot, smile, face) | ✓ | ✓ | ✓ |  |

 EasyChart Polar Plus
 EasyChart EVO
 EasyChart medical LIGHT

| | | | | |
|---|---|---|---|---|
| Dyschromatopsias (numbers) | ✓ | ✓ | ✓ |  |
| Aniseikonia Test | ✓ | ✓ | ✓ |  |
| (in Polar version) | | | |  |
| Crepuscular Visual Acuity | ✓ | ✓ | ✓ |  |
| Contrast sensitivity sinusoidal grating (vertical, oblique) | ✗ | ✓ | ✓ |  |
| Interactive Staircase Albini E | ✗ | ✓ | ✓ |  |
| Interactive Staircase Landolt C | ✗ | ✓ | ✓ |  |
| High crowding Visual Acuity | ✗ | ✓ | ✓ |  |
| Fixation Disparity | ✗ | ✓ | ✓ |  |
| (in Polar version) horizontal, vertical | | | |  |

 EasyChart Polar Plus
 EasyChart EVO
 EasyChart medical LIGHT

| | | | | |
|---|---|---|---|---|
| Reversed contrast Visual Acuity | ✗ | ✓ | ✓ |  |
| Numbers Visual Acuity | ✗ | ✓ | ✓ |  |
| Balancing test (+ Crowen test with letters of different dimensions) | ✗ | ✗ | ✓ |  |
| Retinal Rivalry and Dominance (central and peripheral) | ✗ | ✗ | ✓ |  |
| MKH (cross, needle, arrows, hooks, osteberg, Crowen) | ✗ | ✗ | ✓ |  |
| Stereoscopic Tables (both slashes and rings) | ✗ | ✗ | ✓ |  |
| Dynamic Stereoscopic Tests (rings) | ✗ | ✗ | ✓ |  |



COD. 51-B | 236 SET

Flat sphere positive and negative trial lenses:

0.12 - 0.25 - 0.50 - 0.75 - 1.00 - 1.25 - 1.50 - 1.75 - 2.00 - 2.25
2.50 - 2.75 - 3.00 - 3.25 - 3.50 - 3.75 - 4.00 - 4.50 - 5.00 - 5.50
6.00 - 6.50 - 7.00 - 7.50 - 8.00 - 9.00 - 10.00 - 11.00 - 12.00
13.00 - 14.00 - 15.00 - 16.00 - 18.00 - 20.00

Flat cylinder positive and negative trial lenses:

0.12 - 0.25 - 0.50 - 0.75 - 1.00 - 1.25 - 1.50 - 1.75 - 2.00 - 2.25
2.50 - 2.75 - 3.00 - 3.25 - 3.50 - 4.00 - 4.50 - 5.00 - 6.00

Prism trial lenses:

1DP - 2DP - 3DP - 4DP - 5DP - 6DP - 7DP - 8DP - 9DP - 10DP

Accessory lenses:

Maddox - Cross Line (2 pieces) - Pin Hole - Occluder - Slit -
Red Filter - Green Filter - Plane Lens - Frosted Lens



COD. 51-B | 236 B TRIAL LENSES

in ABS carrying case



COD. 51-BW | 198 B TRIAL LENSES SET

in ABS carrying case w/o positive cylinders

Dimensions: 52,5x41x8,5 cm

COD. 51-BD | 236 B TRIAL LENSES SET

in plastic tray w/o case



COD. 51-BDW | 198 B TRIAL LENSES SET

in plastic tray w/o case and w/o positive cylinders

Dimensions: 51x36x2,5

COD. 51-BM | 236 B TRIAL LENSES SET

in lacquered beech wood case



COD. 51-BMW | 198 B TRIAL LENSES SET

in lacquered beech wood case w/o positive cylinders

Dimensions: 54x39x7 cm

COD. 51-BE | 236 B TRIAL LENSES SET

in sliding cover case




COD. 51-BEW | 198 B TRIAL LENSES SET

in sliding cover case w/o positive cylinders

Dimensions: 56x26x42 cm

COD. 51-BP | 236 B TRIAL LENSES SET

in hanging sliding cover case



COD. 51-BPW | 198 B TRIAL LENSES SET

in hanging sliding cover case w/o positive cylinders

Dimensions: 54x48,5x21 cm

COD. 51-A | 148 SET

Flat sphere positive and negative trial lenses:

0.25 - 0.50 - 0.75 - 1.00 - 1.25 - 1.50 - 1.75 - 2.00 - 2.25 - 2.50 - 2.75 - 3.00 - 3.25
3.50 - 4.00 - 4.50 - 5.00 - 6.00 - 7.00 - 8.00 - 9.00 - 10.00 - 11.00 - 12.00

Flat cylinder positive and negative trial lenses:

0.25 - 0.50 - 0.75 - 1.00 - 1.25 - 1.50 - 1.75 - 2.00 - 2.50 - 3.00

Prism trial lenses:

1DP - 2DP - 3DP - 4DP - 5DP

Accessory lenses:

Maddox - Cross Line - Pin Hole - Occluder - Slit
Red Filter - Green Filter



COD. 51-A | 148 B TRIAL LENSES SET

in ABS carrying case



COD. 51-AW | 128 B TRIAL LENSES SET

in ABS carrying case w/o positive cylinders

Dimensions: 52,5x30x8,5 cm



COD. 51-AD | 148 B TRIAL LENSES SET

in plastic tray w/o case



COD. 51-ADW | 128 B TRIAL LENSES SET

in plastic tray w/o case and w/o positive cylinders

Dimensions: 50,5x25,3x3 cm



COD. 51-AL | 148 B TRIAL LENSES SET

in wood tray



Dimensions: 33x29x3,5 cm





COD. 61-B | 236 SET

Meniscus sphere positive and negative trial lenses:

0.12 - 0.25 - 0.50 - 0.75 - 1.00 - 1.25 - 1.50 - 1.75 - 2.00 - 2.25
2.50 - 2.75 - 3.00 - 3.25 - 3.50 - 3.75 - 4.00 - 4.50 - 5.00 - 5.50
6.00 - 6.50 - 7.00 - 7.50 - 8.00 - 9.00 - 10.00 - 11.00 - 12.00
13.00 - 14.00 - 15.00 - 16.00 - 18.00 - 20.00

Meniscus cylinder positive and negative trial lenses:

0.12 - 0.25 - 0.50 - 0.75 - 1.00 - 1.25 - 1.50 - 1.75 - 2.00 - 2.25
2.50 - 2.75 - 3.00 - 3.25 - 3.50 - 4.00 - 4.50 - 5.00 - 6.00

Prism trial lenses:

1DP - 2DP - 3DP - 4DP - 5DP - 6DP - 7DP - 8DP - 9DP - 10DP

Accessory lenses:

Maddox - Cross Line (2 pieces) - Pin Hole - Occluder - Slit
Red Filter - Green Filter - Plane Lens - Frosted Lens



COD. 61-B | 236 MENISCUS TRIAL LENSES SET

in ABS carrying case



COD. 61-BW | 198 MENISCUS TRIAL LENSES SET

in ABS carrying case w/o positive cylinders

Dimensions: 52,5x41x8,5 cm

COD. 61-BD | 236 MENISCUS TRIAL LENSES SET

in plastic tray w/o case



COD. 61-BDW | 198 MENISCUS TRIAL LENSES SET

in plastic tray w/o case and w/o positive cylinders

Dimensions: 51x36x2,5 cm

COD. 61-BM | 236 MENISCUS TRIAL LENSES SET

in lacquered beech wood case



COD. 61-BMW | 198 MENISCUS TRIAL LENSES SET

in lacquered beech wood case w/o positive cylinders

Dimensions: 54x39x7 cm

COD. 61-BE | 236 MENISCUS TRIAL LENSES SET

in sliding cover case




COD. 61-BEW | 198 MENISCUS TRIAL LENSES SET

in sliding cover case w/o positive cylinders

Dimensions: 56x26x42 cm

COD. 61-BP | 236 MENISCUS TRIAL LENSES SET

in hanging sliding cover case



COD. 61-BPW | 198 MENISCUS TRIAL LENSES SET

in hanging sliding cover case w/o positive cylinders

Dimensions: 54x48,5x21 cm

COD. 61-A | 148 SET

Meniscus sphere positive and negative trial lenses:

0.25 - 0.50 - 0.75 - 1.00 - 1.25 - 1.50 - 1.75 - 2.00 - 2.25 - 2.50
2.75 - 3.00 - 3.25 - 3.50 - 4.00 - 4.50 - 5.00 - 6.00 - 7.00 - 8.00
9.00 - 10.00 - 11.00 - 12.00

Meniscus cylinder positive and negative trial lenses:

0.25 - 0.50 - 0.75 - 1.00 - 1.25 - 1.50 - 1.75 - 2.00 - 2.50 - 3.00

Prism trial lenses:

1DP - 2DP - 3DP - 4DP - 5DP

Accessory lenses:

Maddox - Cross Line - Pin Hole - Occluder
Slit - Red Filter - Green Filter



COD. 61-A | 148 MENISCUS TRIAL LENSES

in ABS carrying case



COD. 61-AW | 128 MENISCUS TRIAL LENSES SET

in ABS carrying case w/o positive cylinders

Dimensions: 52,5x30x8,5 cm



COD. 61-AD | 148 MENISCUS TRIAL LENSES

in plastic tray w/o case



COD. 61-ADW | 128 MENISCUS TRIAL LENSES SET

in plastic tray w/o case and w/o positive cylinders

Dimensions: 50,5x25,3x3 cm




COD.61-AL | 148 MENISCUS TRIAL LENSES

in wood tray

Dimensions: 33x29x3,5 cm





COD. 99-R | SINGLE MENISCUS TRIAL LENS

- meniscus cylinder positive and negative trial lenses:
0.12 - 0.25 - 0.50 - 0.75 - 1.00 - 1.25 - 1.50 - 1.75 - 2.00
2.25 - 2.50 - 2.75 - 3.00 - 3.25 - 3.50 - 4.00 - 4.50 - 5.00 - 6.00
- meniscus sphere positive and negative trial lenses:
0.12 - 0.25 - 0.50 - 0.75 - 1.00 - 1.25 - 1.50 - 1.75 - 2.00 - 2.25
2.50 - 2.75 - 3.00 - 3.25 - 3.50 - 3.75 - 4.00 - 4.50 - 5.00 - 5.50
6.00 - 6.50 - 7.00 - 7.50 - 8.00 - 9.00 - 10.00 - 11.00 - 12.00
13.00 - 14.00 - 15.00 - 16.00 - 18.00 - 20.00






COD. 99-R2 | SINGLE B TRIAL LENS

- flat cylinder positive and negative trial lenses:
0.12 - 0.25 - 0.50 - 0.75 - 1.00 - 1.25 - 1.50 - 1.75 - 2.00 - 2.25
2.50 - 2.75 - 3.00 - 3.25 - 3.50 - 3.75 - 4.00 - 4.50 - 5.00 - 6.00
- flat sphere positive and negative trial lenses:
0.12 - 0.25 - 0.50 - 0.75 - 1.00 - 1.25 - 1.50 - 1.75 - 2.00 - 2.25
2.50 - 2.75 - 3.00 - 3.25 - 3.50 - 3.75 - 4.00 - 4.50 - 5.00
5.50 - 6.00 - 6.50 - 7.00 - 7.50 - 8.00 - 9.00 - 10.00 - 11.00
12.00 - 13.00 - 14.00 - 15.00 - 16.00 - 18.00 - 20.00

COD. 99AC | ACCESSORY LENS

- Maddox - Cross Line - Pin Hole - Occluder - Slit -
Red Filter - Green Filter - Plane Lens - Frosted Lens



COD. 99RP | PRISM TRIAL LENS

- Prism trial lenses 1DP - 2DP - 3DP - 4DP - 5DP - 6DP - 7DP
8DP - 9DP - 10DP





COD. K350 | SINGLE B TRIAL LENS

- flat sphere positive and negative trial lenses:
0.12 - 0.25 - 0.50 - 0.75 - 1.00 - 1.25 - 1.50 - 1.75 - 2.00 - 2.25
2.50 - 2.75 - 3.00 - 3.25 - 3.50 - 4.00 - 4.50 - 5.00 - 6.00
- flat cylinder positive and negative trial lenses:
0.12 - 0.25 - 0.50 - 0.75 - 1.00 - 1.25 - 1.50 - 1.75 - 2.00 - 2.25
2.50 - 2.75 - 3.00 - 3.25 - 3.50 - 3.75 - 4.00 - 4.50 - 5.00
5.50 - 6.00 - 6.50 - 7.00 - 7.50 - 8.00 - 9.00 - 10.00 - 11.00
12.00 - 13.00 - 14.00 - 15.00 - 16.00 - 18.00 - 20.00

COD. K350AC | ACCESSORY LENS

(Maddox - Cross Line - Pin Hole - Occluder - Slit - Red Filter - Green Filter
Plane Lens - Frosted Lens)

COD. K350P | PRISM TRIAL LENS

1DP - 2DP - 3DP - 4DP - 5DP - 6DP - 7DP - 8DP - 9DP - 10DP



COD. 71-C | 236 LENSES SET

In lacquered beech wood case and tray

Flat sphere positive and negative trial lenses:

0.12 - 0.25 - 0.50 - 0.75 - 1.00 - 1.25 - 1.50 - 1.75 - 2.00 - 2.25 - 2.50 - 2.75
 3.00 - 3.25 - 3.50 - 3.75 - 4.00 - 4.50 - 5.00 - 5.50 - 6.00 - 6.50 - 7.00 - 7.50 - 8.00
 9.00 - 10.00 - 11.00 - 12.00 - 13.00 - 14.00 - 15.00 - 16.00 - 18.00 - 20.00

Flat cylinder positive and negative trial lenses:

0.12 - 0.25 - 0.50 - 0.75 - 1.00 - 1.25 - 1.50 - 1.75 - 2.00 - 2.25 - 2.50 - 2.75
 3.00 - 3.25 - 3.50 - 4.00 - 4.50 - 5.00 - 6.00

Prism trial lenses:

1DP - 2DP - 3DP - 4DP - 5DP - 6DP - 7DP - 8DP - 9DP - 10DP

Accessory lenses:

Maddox - Cross Line (2 pieces) - Pin Hole - Occluder - Slit -
 Red Filter - Green Filter - Plane Lens - Frosted Lens

Dimensions: 54x39x7 cm

**COD. 71-P | 68 LENSES SET (PERIMETRY SET)**

In lacquered beech wood case and tray

Flat sphere positive and negative trial lenses:

0.25 - 0.50 - 0.75 - 1.00 - 1.25 - 1.50 - 1.75 - 2.00 - 2.25 - 2.50 - 2.75
 3.00 - 3.50 - 4.00 - 4.50 - 5.00 - 5.50 - 6.00 - 6.50 - 7.00 - 8.00 - 9.00
 10.00 - 11.00 - 12.00 - 13.00 - 14.00 - 16.00

Flat cylinder negative trial lenses:

0.25 - 0.50 - 0.75 - 1.00 - 1.50 - 2.00 - 2.50 - 3.00 - 3.50 - 4.00 - 4.50 - 5.00

Dimensions: 33,5x15,5x6 cm



COD. 78-MT | TABO TRIAL FRAME

Lenses slots: maximum 5 lenses in right slots and 5 lenses in left slots (front ones 3 and back ones 2 for each side)

Interpupillary distance: min. 48 mm to 80 mm




COD. 78-UL | ULTRA LIGHT TRIAL FRAME

Lenses slots: maximum 5 lenses in right slots and 5 lenses in left slots (front ones 3 and back ones 2 for each side)

Binocular interpupillary distance: 50~80 mm

Monocular interpupillary distance: 25~40 mm



COD. RB307-PD56

Trial Frames for Babies aDesigned specifically for pediatric use, ensuring comfort and precise fitting during eye exams.

Size > 41-16 (PD 56)

FEATURES:

- Made of soft and flexible ABS and PP plastic
- No hinges or metal parts
- Small holes at the temple tips
- Adjustable for secure fit
- Elastic band included

INCLUDES:

- Adjustable elastic strap



93-AB | RETINOSCOPY RACKS

In lacquered beech wood case

Green rod: +1D; +2D; +3D; +4D; +5D; +6D; +7D; +8D; +9D; +10D

Slide: +0,5D; +10D

Red rod: -1D; -2D; -3D; -4D; -5D; -6D; -7D; -8D; -9D; -10D

Slide: -0,5D; -10D




93-AF | BERENS PRISM BARS

In lacquered beech wood case


Vertical prism dioptres: 1P; 2P; 3P; 4P; 5P; 6P; 8P; 10P; 12P; 14P; 16P; 18P; 20P; 25P

Horizontal prism dioptres: 1P; 2P; 4P; 6P; 8P; 10P; 12P; 14P; 16P; 18P; 20P; 25P; 30P; 35P; 40P



 **COD. 93-A | FLIPPER TEST**
no lenses and case



 **COD. 93-R | SINGLE LENS FOR RETINOSCOPY RACKS**
positive and negative
0,5D - 1D - 2D - 3D - 4D - 5D - 6D - 7D - 8D - 9D - 10D



COD. 97-LA | SINGLE CROSS CYLINDER

+/- 0.25



COD. 97-MA | SINGLE CROSS CYLINDER

+/- 0.50



COD. 97-LM | CROSS CYLINDERS

+/- 0.25 & +/- 0.50



ASPHERIC LENSES

- Clear and bright image
- Easy-to-hold, high quality lens frame

COD. 20D

Effective diameter: 49 mm
Image magnification: 2.97X
Working distance: 47 mm
Focal length: 51.1 mm
Field of view: 50°
Weight: 53 g



COD. 78D

Effective diameter: 27 mm
Image magnification: 0.77X
Working distance: 10 mm
Focal length: 16.1 mm
Field of view: 74°
Weight: 19 g



COD. 90D

Effective diameter: 20 mm
Image magnification: 0.65X
Working distance: 10 mm
Focal length: 12.5 mm
Field of view: 74°
Weight: 16 g







FE-GROUP
TRUE INNOVATIVE SOLUTIONS

EASYOPHT 
true innovative solutions
division

EasyOpht a brand of FE-Group srl

Via IV Novembre, 118

21058 Solbiate Olona (VA) - Italy

P +39 0331 626490

www.fe-group.it