


# ACCUJECT™




medicel  
swiss evolution

# ACCUJECT™




## The all-rounder

Accuject injectors provides Medical's largest range of applications when injecting IOL's. Its versatile design allows the use of any IOL material and design without limitations on the implant procedure and incision size.



### Accuject (top-loaded / single-use / sterile)

Item Code	Device Name	Incision size (mm)*			IOL type
		wound-assisted	into-the-wound	into-the-bag	
LP604580	3.0-3P Injector Set	2.4	2.8	3.2	3-piece lenses
LP604520	2.5-3P Injector Set	2.2	2.5	2.8-2.9	
LP604570	3.0-1P Injector Set	2.4	2.7-2.8	3.1	1-piece lenses
LP604600	2.4-1P Injector Set	2.1	2.4	2.7	
LP604530	2.2-1P Injector Set	1.9-2.0	2.2	2.5	
LP604590	2.1-1P Injector Set	1.9	2.1-2.2	2.4	
LP604510	2.0-1P Injector Set	1.8-1.9	2.1	2.3-2.4	
LP604540	1.8-1P Injector Set	1.8	2.0	2.2-2.3	
LP604550	1.6-1P Injector Set	1.6	1.8	2.0	hydrophilic only



### Accuject (back-loaded / single-use / sterile)

Item Code	Device Name	Incision size (mm)*			IOL type
		wound-assisted	into-the-wound	into-the-bag	
LP604575	3.0-BL Injector Set	2.4	2.7 - 2.8	3.1	1-piece lenses
LP604605	2.4-BL Injector Set	2.1	2.4	2.7	
LP604535	2.2-BL Injector Set	1.9 - 2.0	2.2	2.5	
LP604595	2.1-BL Injector Set	1.9	2.1 - 2.2	2.4	
LP604515	2.0-BL Injector Set	1.8 - 1.9	2.1	2.3 - 2.4	
LP604545	1.8-BL Injector Set	1.8	2.0	2.2-2.3	

\* Incision sizes are indications only based on average experiences. The individually required incision sizes may vary depending on surgical techniques



instruction for use



medical

Medical AG, Dornierstrasse 11, 9423 Altenrhein, +41 71 727 10 50, medical.com



digital brochure

made in switzerland