

# VeraPlug™ flow

Punctal Occlusion System



The **VeraPlug™ Flow** is designed to provide a simple and effective means of slowing tear drainage.

- Inner channel reduces but does not eliminate tear flow, providing only **partial occlusion**
- Conforms perfectly to the puncta for **unparalleled retention**
- Low-profile dome for **enhanced patient comfort**
- Proprietary shaft design for **easy insertion** and **proper anatomic fit**

## AVAILABLE CONFIGURATIONS

### Sterile Preloaded

Individually packaged and preloaded on a sterile, disposable inserter. 1 pair per box.

SIZE	1 PAIR
Small	VFL-1001
Medium	VFL-1002
Large	VFL-1003
X-Large	VFL-1004

### Nonsterile Bulk

Packaged nonsterile, 1 pair per poly bag, 5 pairs per box. Sterile inserters included.

SIZE	5 PAIRS
Small	BVF-1001
Medium	BVF-1002
Large	BVF-1003
X-Large	BVF-1004

Choose the **sterile preloaded** configuration for the ultimate in convenience or **nonsterile bulk** for increased value.



## SIZE CHART

PUNCTAL OPENING	VERAPLUG™ FLOW SIZE
0.4mm to 0.6mm	Small
0.6mm to 0.7mm	Medium
0.7mm to 0.8mm	Large
0.8mm to 1.0mm	X-Large