

# Classic Trial Lens Set & Frame in an Elegant Carrying Case



## The Reichert Trial Lens Set & Trial Frame

# The Reichert Trial Lens Set & Frame

## Robust and Reliable Lens Set

- Full-view lenses with polished metal rims
- Robust materials and quality craftsmanship for years of reliable performance
- Elegant carrying case



## Comprehensive Set with 228 Lenses

- 70 Spherical Concave Lenses: 0.12 to 20.00 Diopters (0.12 / 0.25D increments)
- 70 Spherical Convex Lenses: 0.12 to 20.00 Diopters (0.12 / 0.25D increments)
- 34 Cylindrical Concave Lenses: 0.12 to 6.00 Diopters (0.12 / 0.25D increments)
- 34 Cylindrical Convex Lenses: 0.12 to 6.00 Diopters (0.12 / 0.25D increments)
- 12 Prismatic Lenses: 0.5(2), 1(2), 2(2), 3, 4, 5, 6, 7, 8
- 8 Accessory Lenses: Red, Green, Slit, Occluder, Pin Hole(2), Cross Cylinder (2)

## Lightweight Trial Frame

A necessary addition to any Trial Lens set, the Reichert Trial Frame provides a comfortable fit and smooth controls for quick, easy lens fitting.



## Familiar Design and Operation

The robust design, high-quality materials, and precision assembly ensure years of trouble-free use.

### Catalog Numbers

15010	TF50 Trial Frame
15015	TS75 Trial Lens Set

*Specifications subject to change without notice*

Call your authorized Reichert Ophthalmic Instruments Distributor today.

Distributor list available at [www.reichertoi.com](http://www.reichertoi.com)

#### Important Notice

For your own protection, only purchase your Reichert Ophthalmic Instruments from, and have them serviced by, an Authorized Reichert Distributor.



Reichert products are designed and manufactured under quality processes meeting ISO 9001 requirements



Reichert, Inc.  
3374 Walden Avenue, Depew, NY 14043  
Toll Free: 888-849-8955 | Tel: (716) 686-4500  
Fax: (716) 686-4545 | [www.reichert.com](http://www.reichert.com)