

PSL one

Excellence  
in optics



Keeler

## Keeler – British manufacturing experience and quality delivering great value

### Precision machined aluminium chassis

Keeler designed and manufactured with the strength to withstand the daily traumas of the demanding ophthalmic environment. Lightweight, durable and robust.

### Advanced optics, 10x magnification

Top quality multi layered glass is crafted into our Keeler multi-element and multi-layer coated optical system to give maximum performance. Clarity is guaranteed.

### Controllable illumination

Proven halogen illumination with illumination levels controlled by a rheostat, giving you total control of your light levels. From zero to maximum, you choose.

### Big Slit Lamp features, Portable usability

Unique features such as the fixation targets, and the 1mm square light patch for assessing anterior chamber flare sets the new standard for Portable Slit Lamps.

### For the digital age

Simply connect an iPhone 4 and the PSL becomes a portable, digital slit lamp. What could be simpler or more cost efficient?



#### Slit wheel



#### Filter wheel

Red free, Blue, Neutral density (0.8) and Clear



#### Magnification



## Slit and filter system

Slit type	Rotating Slit wheel selection Slit Lamp
Slit length	12mm
Slit widths	0.15mm, 0.5mm, 0.8mm and 1.6mm slits, 12mm circle and a 1mm square
Filters	Red free, Blue, Neutral density 0.8 and Clear
IR protection	Inbuilt IR cut filter
Slit angle	+/- 60°
Illumination control	Continuously variable from low to full brightness
Fixation targets	LED system, double click of trigger

### Power supply

Power supply unit	Switch mode, (110V – 240V) +/- 10% multi-plug
-------------------	--

Power supply output	30VA (12V DC 2.5A)
---------------------	--------------------

### Biomicroscope

Type	Binocular Hand held Biomicroscope Slit Lamp
Optics	Converging binoculars at 13°
Magnification	10x
Objective lens working distance @ 10x	16mm

### Ordering details

PSL One	3010-P-2001
iPhone 4 adapter	3010-P-7010 (included)
Carrying case	3010-P-7000

iPhone is a registered trade mark of Apple Inc  
iPhone not included with adapter or PSL

**PSL one**

**Keeler Limited, Clewer Hill Road, Windsor, Berkshire SL4 4AA**

+44 (0)1753 857177 | [www.keeler.co.uk](http://www.keeler.co.uk) | [info@keeler.co.uk](mailto:info@keeler.co.uk)

**Keeler Instruments Inc. 456 Parkway Broomall, PA 19008 USA**

Toll Free: 1 800 523 5620 | Fax: 1 610 353 7814 | [www.keelerusa.com](http://www.keelerusa.com)